

***Product Specifications***

<b>Amount</b>	◆ 1 mg
<b>Form</b>	◆ Lyophilized
<b>Purity</b>	◆ > 90% by SDS-PAGE
<b>Buffer</b>	◆ PBS, pH 7.4 with 5% trehalose
<b>Endotoxin</b>	◆ < 0.10 EU per 1 µg of the antibody by the LAL method.
<b>Preservative</b>	◆ None
<b>Reconstitution</b>	◆ Reconstitute with sterile PBS. If 0.5 mL of PBS is used, the antibody concentration will be 2 mg/mL.
<b>Storage</b>	◆ 2° - 8° C
<b>Source</b>	◆ Normal rat serum