

DESCRIPTION

Source Chinese Hamster Ovary cell line, CHO-derived
Met1-Ile259
Accession # NP_055236

N-terminal Sequence Analysis Ser26

Predicted Molecular Mass 25.8 kDa

SPECIFICATIONS

SDS-PAGE 32-36 kDa, reducing conditions

Activity Measured by its ability to inhibit Wnt induced TCF reporter activity in HEK293 human embryonic kidney cells. Recombinant Human Dkk-2 (Catalog # 6628-DK/CF) inhibits a constant dose of 500 ng/mL of Recombinant Human Wnt-3a (Catalog # 5036-WN). The ED₅₀ for this effect is 50-300 ng/mL.

Endotoxin Level <0.10 EU per 1 µg of the protein by the LAL method.

Purity >95%, by SDS-PAGE visualized with Silver Staining and quantitative densitometry by Coomassie® Blue Staining.

Formulation Lyophilized from a 0.2 µm filtered solution in PBS and EDTA. See Certificate of Analysis for details.

PREPARATION AND STORAGE


Reconstitution Reconstitute at 100 µg/mL in PBS.

Shipping The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

Stability & Storage Use a manual defrost freezer and avoid repeated freeze-thaw cycles.

- 12 months from date of receipt, -20 to -70 °C as supplied.
- 1 month, 2 to 8 °C under sterile conditions after reconstitution.
- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

DATA



BACKGROUND

Dickkopf related protein 2 (Dkk-2) is a member of the Dickkopf family of secreted Wnt modulators (1-3). Dkk proteins contain a signal peptide and two conserved cysteine-rich domains that are separated by a linker region. The second cysteine-rich domain mediates Dkk-2 binding activities, and its interaction with β -propeller domains of LRP-5/6 has been mapped (2-4, 7). The 226 amino acid (aa), ~35 kDa mature human Dkk-2 shares 96%, 97%, 97%, 97%, 97% and 98% aa identity with mouse, rat, canine, equine, bovine and porcine Dkk-2, respectively. Mouse Dkk-2 can activate the canonical Wnt signaling pathway in *Xenopus* embryos, showing evolutionary conservation of function (5). Dkk proteins modify Wnt engagement of a receptor complex composed of a Frizzled protein and a low-density lipoprotein receptor-related protein, either LRP-5 or LRP-6 (3). Also, Kremen-1 and Kremen-2 are high affinity receptors for Dkk-1 and Dkk-2 (9). When LRP-6 is over-expressed, direct high-affinity binding of Dkk-2 to LRP can enhance canonical Wnt signaling (6-8). However, when Dkk-2 and LRP-6 form a ternary complex with Kremen-2, Wnt signaling is inhibited due to internalization of Dkk-2/LRP6/Krm2 complexes (9, 10). Thus, depending on the cellular context, Dkk-2 can either activate or inhibit canonical Wnt signaling (3). In contrast, binding of Dkk-1 or Dkk-4 to LRP is consistently antagonistic (3). Dkk proteins are expressed in mesenchymal tissues and control epithelial transformations. Dkk-2 expression has been studied most in bone and eye, although it is expressed as early as periimplantation in mice (11). Mouse Dkk-1 or Dkk-2 deficiencies have opposite effects on bone homeostasis, despite down-regulating Wnt antagonism in both cases (12, 13). Dkk-2 expression is induced by Wnts in bone, and is thought to enhance bone density by promoting terminal differentiation of osteoblasts and mineral deposition (12). In contrast, Dkk-1 negatively regulates late osteoblast proliferation, which limits bone density (13). Dkk-2-deficient mice are blind, exhibiting faulty differentiation of corneal epithelium and ectopic blood vessels in the periocular mesenchyme (14, 15).

References:

1. Monaghan, A.P. *et al.* (1999) *Mech. Dev.* **87**:45.
2. Krupnik, V.E. *et al.* (1999) *Gene* **238**:301.
3. Niehrs, C. (2006) *Oncogene* **25**:7469.
4. Chen, L. *et al.* (2008) *J. Biol. Chem.* **283**:23364.
5. Wu, W. *et al.* (2000) *Current Biol.* **10**:1611.
6. Mao, B. *et al.* (2001) *Nature* **411**:321.
7. Li, L. *et al.* (2002) *J. Biol. Chem.* **277**:5977.
8. Brott, B. and S.Y. Sokol (2002) *Mol. Cell. Biol.* **22**:6100.
9. Mao, B. *et al.* (2002) *Nature* **417**:664.
10. Mao, B. and C. Niehrs (2003) *Gene* **302**:179.
11. Zhang, Y. *et al.* (2009) *J. Reprod. Dev.* **55**:17.
12. Li, X. *et al.* (2005) *Nat. Genet.* **37**:945.
13. van der Horst, G. *et al.* (2005) *J. Bone Miner. Res.* **20**:1867.
14. Mukhopadhyay, M. *et al.* (2006) *Development* **133**:2149.
15. Gage, P.J. *et al.* (2008) *Dev. Biol.* **317**:310.