

DESCRIPTION

Source *E. coli*-derived
Ala21-Ser153, with an N-terminal Met
Accession # NP_001156652

N-terminal Sequence Analysis Met & Ala21

Structure / Form Monomer

Predicted Molecular Mass 15.2 kDa

SPECIFICATIONS

SDS-PAGE 14 kDa, reducing conditions

Activity Measured in a cell proliferation assay using CTLL-2 mouse cytotoxic T cells.
The ED₅₀ for this effect is 0.015-0.075 ng/mL.

Endotoxin Level <0.01 EU per 1 µg of the protein by the LAL method.

Purity >95%, by SDS-PAGE visualized with Silver Staining and quantitative densitometry by Coomassie® Blue Staining.

Formulation Lyophilized from a 0.2 µm filtered solution in Sodium Acetate with BSA as a carrier protein. See Certificate of Analysis for details.

PREPARATION AND STORAGE

Reconstitution Reconstitute at 100 µg/mL in 100 mM acetic acid containing at least 0.1% human or bovine serum albumin.

Shipping The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.

Stability & Storage Use a manual defrost freezer and avoid repeated freeze-thaw cycles.

- 12 months from date of receipt, -20 to -70 °C as supplied.
- 1 month, 2 to 8 °C under sterile conditions after reconstitution.
- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

BACKGROUND

Interleukin-2 (IL-2) is a 15 kDa, O-glycosylated, four α-helix bundle cytokine that has potent stimulatory activity for antigen-activated T cells. It is expressed by CD4⁺ and CD8⁺ T cells, γδ T cells, B cells, dendritic cells, and eosinophils (1-3). Mature rabbit IL-2 shares 79%, 54%, and 63% aa sequence identity with human, mouse, and rat IL-2, respectively. Human and mouse IL-2 exhibit cross-species activity (4). The receptor for IL-2 consists of three subunits that are present on the cell surface in varying preformed complexes (5-7). The 55 kDa IL-2 Rα is specific for IL-2 and binds with low affinity. The 75 kDa IL-2 Rβ, which is also a component of the IL-15 receptor, binds IL-2 with intermediate affinity. The 64 kDa common gamma chain γc/IL-2 Rγ, which is shared with the receptors for IL-4, -7, -9, -15, and -21, does not independently interact with IL-2. Upon ligand binding, signal transduction is performed by both IL-2 Rβ and γc. IL-2 is best known for its autocrine and paracrine activity on T cells. It drives resting T cells to proliferate and induces IL-2 and IL-2 Rα synthesis (1, 2). It contributes to T cell homeostasis by promoting the Fas-induced death of naïve CD4⁺ T cells but not activated CD4⁺ memory lymphocytes (8). IL-2 plays a central role in the expansion and maintenance of regulatory T cells, although it inhibits the development of Th17 polarized cells (9-11). Thus, IL-2 may be a key cytokine in the natural suppression of autoimmunity (12-14).

References:

1. Ma, A. *et al.* (2006) *Annu. Rev. Immunol.* **24**:657.
2. Gaffen, S.L. and K.D. Liu (2004) *Cytokine* **28**:109.
3. Perkins, H.D. *et al.* (2000) *Cytokine* **12**:555.
4. Mosmann, T.R. *et al.* (1987) *J. Immunol.* **138**:1813.
5. Liparoto, S.F. *et al.* (2002) *Biochemistry* **41**:2543.
6. Wang, X. *et al.* (2005) *Science* **310**:1159.
7. Bodnar, A. *et al.* (2008) *Immunol. Lett.* **116**:117.
8. Jaleco, S. *et al.* (2003) *J. Immunol.* **171**:61.
9. Malek, T.R. (2003) *J. Leukoc. Biol.* **74**:961.
10. Laurence, A. *et al.* (2007) *Immunity* **26**:371.
11. Kryczek, I. *et al.* (2007) *J. Immunol.* **178**:6730.
12. Afzali, B. *et al.* (2007) *Clin. Exp. Immunol.* **148**:32.
13. Dooms, H. and A.K. Abbas (2010) *Eur. J. Immunol.* **40**:1538.
14. Liu, R. *et al.* (2010) *Eur. J. Immunol.* **40**:1577.