

DESCRIPTION


Specificity	Detects TGF- β 2/1.2 in direct ELISAs and Western blots. In direct ELISAs and Western blots (non-reducing conditions), less than 2% cross-reactivity with TGF- β 1, TGF- β 3, and TGF- β 5 is observed.
Source	Polyclonal Goat IgG
Purification	Antigen Affinity-purified
Immunogen	Porcine platelet-derived TGF- β 2
Endotoxin Level	<0.10 EU per 1 μ g of the antibody by the LAL method.
Formulation	Lyophilized from a 0.2 μ m filtered solution in PBS with Trehalose. See Certificate of Analysis for details. *Small pack size (-SP) is supplied as a 0.2 μ m filtered solution in PBS.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

	Recommended Concentration	Sample
Western Blot	0.1 μ g/mL	Recombinant Human TGF- β 2 (Catalog # 302-B2) Recombinant Human TGF- β 1.2 (Catalog # 304-B3)
Immunohistochemistry	5-15 μ g/mL	Immersion fixed paraffin-embedded sections of human prostate cancer tissue
Neutralization	Measured by its ability to neutralize TGF- β 2 inhibition of IL-4-dependent proliferation in the HT-2 mouse T cell line. Tsang, M. <i>et al.</i> (1995) Cytokine 7:389. The Neutralization Dose (ND ₅₀) is typically 0.01-0.03 μ g/mL in the presence of 1 ng/mL Porcine TGF- β 2 and 7.5 ng/mL Recombinant Mouse IL-4.	

DATA



PREPARATION AND STORAGE

Reconstitution	Reconstitute at 0.2 mg/mL in sterile PBS.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below. *Small pack size (-SP) is shipped with polar packs. Upon receipt, store it immediately at -20 to -70 °C
Stability & Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <ul style="list-style-type: none"> • 12 months from date of receipt, -20 to -70 °C as supplied. • 1 month, 2 to 8 °C under sterile conditions after reconstitution. • 6 months, -20 to -70 °C under sterile conditions after reconstitution.