

DESCRIPTION

Species Reactivity	Human
Specificity	Detects human IL1RAPL2 in direct ELISAs and Western blots. In direct ELISAs, approximately 35% cross-reactivity with recombinant mouse IL1RAPL2 is observed. In Western blots, approximately 5% cross-reactivity with recombinant human (rh) IL-1&nbs
Source	Polyclonal Goat IgG
Purification	Antigen Affinity-purified
Immunogen	Mouse myeloma cell line NS0-derived recombinant human IL1RAPL2 Thr17-Glu356 Accession # Q9NP60
Conjugate	Alexa Fluor 488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

CyTOF-ready	Optimal dilution of this antibody should be experimentally determined.
Western Blot	Optimal dilution of this antibody should be experimentally determined.
Flow Cytometry	Optimal dilution of this antibody should be experimentally determined.
Immunohistochemistry	Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

The Interleukin 1 receptor family (IL-1 R) comprises at least eleven members including IL-1 RI (IL-1 R1), IL-1 RII (IL-1 R2), IL-1 RAcP (IL-1 R3), ST2 (T1/IL-1 R4), IL-18 Rα (IL-1 Rrp/IL-1 R5), IL-1 Rrp2 (IL-1 RL2/IL-1 R6), IL-18 Rβ (AcPL/IL-1 R7), IL1RAPL1 (TIGIRR-2/IL-1 R8), and IL1RAPL2 (TIGIRR-1/IL-1 R9) (1). All family members possess three immunoglobulin (Ig)-like domains in their extracellular region. Most members also have an intracellular TIR (Toll-like receptor/IL-1 receptor signaling) domain that is also conserved in the Toll-like receptor family. Related proteins, SIGIRR (single Ig domain-containing IL-1 R-related molecule) and IL-18BP, differ from the other members by having only one Ig domain (1). IL-1 receptor accessory protein-like 2 (IL-1 RAPL2) is alternately known as IL-1 R9 and three immunoglobulin domain containing IL-1 receptor-related molecule 1 (TIGIRR-1) and is expressed in the brain (2). Its sequence predicts an 686 amino acid (aa) residue type I transmembrane glycoprotein with a 17 aa signal peptide, a 339 aa extracellular region containing three Ig-like domains, an 18 aa transmembrane domain and a 312 aa cytoplasmic tail (3). By comparison to other IL-1 receptor family proteins, IL1RAPL2 has a C-terminal cytoplasmic extension beyond the TIR domain that is found in IL1RAPL1 and SIGIRR but not other family members (3). Human and mouse IL1RAPL2 share approximately 95% aa sequence identity. Human IL1RAPL2 is most homologous (63%) to IL1RAPL1, a receptor protein that is highly expressed in hippocampus and is involved in X-linked mental retardation (4, 5). Genes for both have been localized to human chromosome Xq22. A ligand for IL1RAPL2 has not been identified (1).

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.