

DESCRIPTION

Species Reactivity	Human
Specificity	Detects human TWEAK in direct ELISAs and Western blots. In direct ELISAs, less than 5% cross-reactivity with recombinant human (rh) BAFF, rhFas Ligand, rhGITR Ligand, rhOX40 Ligand, rhTRAIL and less than 1% cross-reactivity with rhAPRIL, rhLIGHT, rhT
Source	Polyclonal Goat IgG
Purification	Antigen Affinity-purified
Immunogen	<i>E. coli</i> -derived recombinant human TWEAK Arg93-His249 Accession # Q4ACW9
Conjugate	Alexa Fluor 700 Excitation Wavelength: 675-700 nm Emission Wavelength: 723 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

Neutralization	Optimal dilution of this antibody should be experimentally determined.
Western Blot	Optimal dilution of this antibody should be experimentally determined.
Immunocytochemistry	Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

TWEAK is a type II transmembrane protein belonging to the TNF superfamily (1). It contains a short cytoplasmic domain (amino acids (aa) 1-18), the transmembrane domain (aa 19-42) and an extracellular domain (aa 43-249). The extracellular domains of human and murine TWEAK share 89% aa sequence identity. A soluble form of TWEAK is generated from the membrane-associated molecules by proteolytic cleavage after Arg 93 suggesting that TWEAK may have long-range effects. TWEAK is expressed widely in many tissues and cells. At least two receptors that bind TWEAK have been identified (2-4). Death Receptor 3 (DR3), also known as TNFRSF12, Apo-3, LARD, WSL-1, or TRAMP, is a TNF receptor superfamily member that is expressed predominantly in tissues with high lymphocyte content (2). It has been suggested that induction of cell death by TWEAK-DR3 interaction involves the activation of NF-κB. In cells that lack DR3, alternate pathways of TWEAK-induced cell death mediated by receptors distinct from DR3 have been suggested (5, 6). TWEAK receptor (TWEAKR, alternatively known as FN14), is a novel TNF receptor superfamily member that also binds TWEAK (3, 4). It is a mitogen-inducible gene that is expressed in fibroblasts, hepatocellular carcinomas and endothelial cells. TWEAK-TWEAKR interaction has been shown to play a role in endothelial cell growth and migration. This effect of TWEAK is not mediated by an up-regulation of VEGF (7).

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.