

#### DESCRIPTION

|                                                                                                                                                                                                         |                                                                                                                                                                        |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Species Reactivity</b>                                                                                                                                                                               | Human                                                                                                                                                                  |
| <b>Specificity</b>                                                                                                                                                                                      | Detects human Semaphorin 6A in direct ELISAs and Western blots. In direct ELISAs, approximately 40% cross-reactivity with recombinant mouse Semaphorin 6A is observed. |
| <b>Source</b>                                                                                                                                                                                           | Polyclonal Goat IgG                                                                                                                                                    |
| <b>Purification</b>                                                                                                                                                                                     | Antigen Affinity-purified                                                                                                                                              |
| <b>Immunogen</b>                                                                                                                                                                                        | Mouse myeloma cell line NS0-derived recombinant human Semaphorin 6A<br>Gly19-Thr649<br>Accession # Q9H2E6                                                              |
| <b>Conjugate</b>                                                                                                                                                                                        | Alexa Fluor 532<br>Excitation Wavelength: 534 nm<br>Emission Wavelength: 553 nm                                                                                        |
| <b>Formulation</b>                                                                                                                                                                                      | Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide                                                                                                           |
| *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions. |                                                                                                                                                                        |

#### APPLICATIONS

**Please Note:** Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

|                             |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| <b>CyTOF-ready</b>          | Optimal dilution of this antibody should be experimentally determined. |
| <b>Western Blot</b>         | Optimal dilution of this antibody should be experimentally determined. |
| <b>Flow Cytometry</b>       | Optimal dilution of this antibody should be experimentally determined. |
| <b>Immunohistochemistry</b> | Optimal dilution of this antibody should be experimentally determined. |

#### PREPARATION AND STORAGE

|                                |                                                                                                                   |
|--------------------------------|-------------------------------------------------------------------------------------------------------------------|
| <b>Shipping</b>                | The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below. |
| <b>Stability &amp; Storage</b> | Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied                          |

#### BACKGROUND

The semaphorins constitute a large family of secreted, glycosylphosphatidylinositol (GPI)-anchored and transmembrane cell signaling molecules. Depending on their domain organization and species origin, these proteins can be classified into eight groups. To date, at least 19 vertebrate Semaphorins belonging to five groups (class 3 to 7), have been identified. All Semaphorins contain a conserved 500 amino acid (aa) Sema domain at the amino-terminus. Semaphorins are best known for their roles in axon guidance during neuronal development. They are also expressed in non-neuronal tissues and are involved in angiogenesis, hematopoiesis, organogenesis, and the regulation of immune functions (1, 2). Class 6 Semaphorins (Sema 6) are transmembrane proteins that share homology with the axon-guiding insect Sema 1A. Human Sema 6A (Vla) cDNA predicts a 1,030 aa protein comprised of an extracellular domain, a transmembrane domain, and a long cytoplasmic tail (3, 4). A secreted form of Sema 6A can repel sympathetic and dorsal root ganglion axons *in vitro*, indicating a traditional role as an axon guidance signal (5). There is evidence, however, that Sema 6A also functions as a guidance receptor. Sema 6A mutants show a defect in thalamocortical neuron projection that is cell autonomous, and the cytoplasmic tails for Sema 6 contain binding sites for intracellular regulatory molecules such as Evl and Src (6).

#### PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.