

**DESCRIPTION**

<b>Species Reactivity</b>	Human
<b>Specificity</b>	Detects human PACT in direct ELISAs and Western blots.
<b>Source</b>	Polyclonal Goat IgG
<b>Purification</b>	Antigen Affinity-purified
<b>Immunogen</b>	<i>E. coli</i> -derived recombinant human PACT Ser2-Lys313 Accession # O75569
<b>Conjugate</b>	Alexa Fluor Plus 555 Excitation Wavelength: 558 nm Emission Wavelength: 572 nm
<b>Formulation</b>	Supplied 0.2 mg/mL in a saline solution containing BSA and Sodium Azide.  *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

**APPLICATIONS**

**Please Note:** Optimal dilutions should be determined by each laboratory for each application. General Protocols are available in the Technical Information section on our website.

**Western Blot** Optimal dilution of this antibody should be experimentally determined.

**DATA**

**PREPARATION AND STORAGE**

**Shipping** The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.

**Stability & Storage** Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

**BACKGROUND**

PACT (PKR Activator; also HSD14, PRKRA and RAX [in mouse]) is a 34-37 kDa member of the PRKRA family of molecules. It is a ubiquitously expressed, cytoplasmic, 313 amino acid (aa) protein that contains three discrete double-stranded (ds) RNA binding motifs (aa 31-101; 126-194; 240-308). It exists as a monomer, homodimer, and heterodimer with at least two partners, TRBP and Dicer. Two isoform variants are known to exist, including one that utilizes an alternative start site at Met26, and another that shows a 10 aa substitution for aa 1-21. PACT is perhaps best known for its relationship with PKR, a 52-58 kDa Ser/Thr kinase that downmodulates protein synthesis and initiates apoptosis. PKR is constitutively silent, and activated via two principal ways; the first is through binding to (viral) dsRNA, and the second is through binding to an activator, PACT. PACT is kept in check upstream through dimerization with at 48-49 kDa regulatory protein termed TRBP. Following growth factor withdrawal or exposure to bacterial products such as LPS, PACT is phosphorylated at Ser287, and this inactive PACT:TRBP heterodimer dissociates, freeing up PACT to bind to PKR. PACT's interaction with PKR induces an autophosphorylation event on PKR. This allows for an activated PKR interaction with eIF2a, causing its phosphorylation and a subsequent reduction of protein synthesis. PACT may also act independently of PKR, and serve as a component of RISC, a complex that downregulates mRNA translation. Within this complex, PACT likely associates with Dicer in concert with Argonaute. Full-length PACT shares 98% aa sequence identity with the mouse ortholog to human PACT termed RAX.

**PRODUCT SPECIFIC NOTICES**

This product is provided under an intellectual property license from Life Technologies Corporation. The transfer of this product is conditioned on the buyer using the purchased product solely in research conducted by the buyer, excluding contract research or any fee for service research, and the buyer must not (1) use this product or its components for (a) diagnostic, therapeutic or prophylactic purposes; (b) testing, analysis or screening services, or information in return for compensation on a per-test basis; or (c) manufacturing or quality assurance or quality control, and/or (2) sell or transfer this product or its components for resale, whether or not resold for use in research. For information on purchasing a license to this product for purposes other than as described above, contact Life Technologies Corporation, 5781 Van Allen Way, Carlsbad, CA 92008 USA or [outlicensing@thermofisher.com](mailto:outlicensing@thermofisher.com).