

## DESCRIPTION

|                                                                                                                                                                                                         |                                                                                     |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| <b>Species Reactivity</b>                                                                                                                                                                               | Human                                                                               |
| <b>Specificity</b>                                                                                                                                                                                      | Detects human PACT in direct ELISAs and Western blots.                              |
| <b>Source</b>                                                                                                                                                                                           | Polyclonal Goat IgG                                                                 |
| <b>Purification</b>                                                                                                                                                                                     | Antigen Affinity-purified                                                           |
| <b>Immunogen</b>                                                                                                                                                                                        | <i>E. coli</i> -derived recombinant human PACT<br>Ser2-Lys313<br>Accession # O75569 |
| <b>Conjugate</b>                                                                                                                                                                                        | Alexa Fluor 647<br>Excitation Wavelength: 650 nm<br>Emission Wavelength: 668 nm     |
| <b>Formulation</b>                                                                                                                                                                                      | Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide                        |
| *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions. |                                                                                     |

## APPLICATIONS

**Please Note:** Optimal dilutions should be determined by each laboratory for each application. General Protocols are available in the Technical Information section on our website.

**Western Blot** Optimal dilution of this antibody should be experimentally determined.

## PREPARATION AND STORAGE

|                                |                                                                                                                   |
|--------------------------------|-------------------------------------------------------------------------------------------------------------------|
| <b>Shipping</b>                | The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below. |
| <b>Stability &amp; Storage</b> | Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied                          |

## BACKGROUND

PACT (PKR Activator; also HSD14, PRKRA and RAX [in mouse]) is a 34-37 kDa member of the PRKRA family of molecules. It is a ubiquitously expressed, cytoplasmic, 313 amino acid (aa) protein that contains three discrete double-stranded (ds) RNA binding motifs (aa 31-101; 126-194; 240-308). It exists as a monomer, homodimer, and heterodimer with at least two partners, TRBP and Dicer. Two isoform variants are known to exist, including one that utilizes an alternative start site at Met26, and another that shows a 10 aa substitution for aa 1-21. PACT is perhaps best known for its relationship with PKR, a 52-58 kDa Ser/Thr kinase that downmodulates protein synthesis and initiates apoptosis. PKR is constitutively silent, and activated via two principal ways; the first is through binding to (viral) dsRNA, and the second is through binding to an activator, PACT. PACT is kept in check upstream through dimerization with at 48-49 kDa regulatory protein termed TRBP. Following growth factor withdrawal or exposure to bacterial products such as LPS, PACT is phosphorylated at Ser287, and this inactive PACT:TRBP heterodimer dissociates, freeing up PACT to bind to PKR. PACT's interaction with PKR induces an autophosphorylation event on PKR. This allows for an activated PKR interaction with eIF2α, causing its phosphorylation and a subsequent reduction of protein synthesis. PACT may also act independently of PKR, and serve as a component of RISC, a complex that downregulates mRNA translation. Within this complex, PACT likely associates with Dicer in concert with Argonaute. Full-length PACT shares 98% aa sequence identity with the mouse ortholog to human PACT termed RAX.

## PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc., and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.