

DESCRIPTION

Species Reactivity	Mouse
Specificity	Detects mouse JMJD6/PSR in direct ELISAs and Western blots.
Source	Polyclonal Goat IgG
Purification	Antigen Affinity-purified
Immunogen	Mouse myeloma cell line NS0-derived recombinant mouse JMJD6/PSR Ala289-Arg403 Accession # Q9ERI5
Conjugate	Alexa Fluor 750 Excitation Wavelength: 749 nm Emission Wavelength: 775 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide
*Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.	

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. General Protocols are available in the Technical Information section on our website.

Western Blot Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

JMJD6 (Jumonji [Japanese for "cruciform/cross"] domain-containing protein 6; also PTDSR/PSR and Lysyl-hydroxylase JMJD6) is a 47-55 kDa member of the JMJD6 family, cupin superfamily of proteins. It is induced by hypoxia, and expressed in skeletal muscle, macrophages and possibly endothelial cells. JMJD6 is reported to demethylate histones (which it does not), serve as a phosphatidyl serine receptor (which it might), and act as a lysyl hydroxylase (which it does). In the nucleus, it appears to regulate pre-RNA splicing by hydroxylating various splicing factors such as U2AF65. Mouse JMJD6 is 403 amino acids (aa) in length. It contains five NLS's plus a JmjC domain (aa 141-305) and three iron-based catalytic sites. There are at least three potential isoform variants. One contains a 29 aa substitution for aa 270-403, a second possesses a 30 aa substitution for aa 302-314, and a third shows a 21 aa substitution for aa 315-403. JMJD6 is known to form high MW homo-trimers, -tetramers and -oligomers. Full-length mouse JMJD6 shares 98% and 99% aa sequence identity with human and rat JMJD6, respectively.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.