

DESCRIPTION

Species Reactivity	Human/Mouse
Specificity	Detects human NeuroD2 in direct ELISAs and Western blots. In direct ELISAs, approximately 20% cross-reactivity with recombinant human NeuroD1 is observed.
Source	Polyclonal Goat IgG
Purification	Antigen Affinity-purified
Immunogen	<i>E. coli</i> -derived recombinant human NeuroD2 Ser174-Asn382 Accession # Q15784
Conjugate	Alexa Fluor 532 Excitation Wavelength: 534 nm Emission Wavelength: 553 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

Western Blot Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.

Stability & Storage Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

NeuroD2 (Neurogenic differentiation factor 2; also NDRF and KW8) is a 41 kDa nuclear member of the neuroD family of transcription factors. It is expressed in developing and mature neurons such as hippocampal granular neurons, and acts (in part) to repress factors that would otherwise block multipotential cell commitment to a neuronal lineage. NeuroD2 is presumed to act as a heterodimer with other bHLH transcription factors. Human NeuroD2 is 382 amino acids (aa) in length. It contains a poly-Glu region (aa 82-91), an NLS (aa 107-113), a DNA-binding HLH domain (aa 119-178) and a poly-Gly segment (aa 282-285). Based on rat studies, human NeuroD2 will undergo variable phosphorylation in the C-terminal region. Over aa 174-382, human and mouse NeuroD2 possess identical amino acid sequences.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.