

DESCRIPTION

Species Reactivity	Human
Specificity	Detects human TRAF-5 in direct ELISAs. In direct ELISAs, approximately 15% cross-reactivity with recombinant human (rh) TRAF-3 is observed, and less than 5% cross-reactivity with rhTRAF-2 and rhTRAF-4 is observed.
Source	Polyclonal Sheep IgG
Purification	Antigen Affinity-purified
Immunogen	<i>E. coli</i> -derived recombinant human TRAF-5 Lys353-Leu557 Accession # O00463
Conjugate	Alexa Fluor 488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

Immunohistochemistry Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.

Stability & Storage Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

TRAF-5 (TNF Receptor-Associated Factor 5; also RNF84) is a 60-66 kDa member of the A subfamily of the TRAF family of molecules. It is widely expressed, and serves as an adaptor molecule that links select TNFR family members with downstream signaling elements. In particular, it interacts with LTbR, TNFR2, OX40, CD27, CD40, CD30, RANK and TAJ/TROY, and induces NFκB activation. It forms homotrimers and heteromers with TRAF3. Human TRAF5 is 557 amino acids (aa) in length. It contains one RING-type Zn finger region (aa 45-85), two TRAF-type Zn finger regions (aa 127-181 and 182-239), a coiled-coil region (aa 237-342) and a C-terminal protein-protein interaction MATH domain (aa 403-549). There is one potential alternative start site 64 aa upstream of the standard start site. Full-length human TRAF-5 shares 80% aa sequence identity with mouse TRAF-5.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.