

DESCRIPTION

Species Reactivity	Human
Specificity	Detects human SODD/BAG4 in Western blots.
Source	Polyclonal Goat IgG
Purification	Antigen Affinity-purified
Immunogen	<i>E. coli</i> -derived recombinant human SODD/BAG4 Glu122-Leu457 Accession # O95429
Conjugate	Alexa Fluor 532 Excitation Wavelength: 534 nm Emission Wavelength: 553 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide
*Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.	

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

Knockout Validated	Optimal dilution of this antibody should be experimentally determined.
Western Blot	Optimal dilution of this antibody should be experimentally determined.
Immunocytochemistry	Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

Tumor necrosis factor receptor-1 (TNF-R1) and other TNF receptor superfamily members, such as DR3, contain intracellular death domains (DD) and are capable of initiating apoptosis when activated by their ligands. Silencer of Death Domains (SODD) was identified as being involved in the cellular mechanism to protect against ligand-independent signaling by TNF-R1 and other DD receptors. SODD, also known as Bcl-2-Associated Athanogene 4 (BAG4), is a 457 amino acid (aa), anti-apoptotic protein that functions through interactions with a variety of proteins including Bcl-2, Raf-protein kinase, steroid hormone receptors, growth factor receptors, and members of the heat shock protein 70 kDa family. SODD is a ubiquitously expressed, cytoplasmic protein that contains a C-terminal BAG domain that can bind and inhibit the chaperone activity of Hsc70/Hsp70. The association of SODD with the DD of TNF-R1 prevents constitutive activation of the TNF-R1 signaling pathway. Binding of TNF to TNF-R1 releases SODD and permits adapter molecules such as TRADD to associate with TNF-R1 leading to the activation of TNF signaling pathways such as apoptosis and NFκB activation.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.