

## DESCRIPTION

|                           |                                                                                                                                                                                                                                                                             |
|---------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Species Reactivity</b> | Human                                                                                                                                                                                                                                                                       |
| <b>Specificity</b>        | Detects human MDGA1 in direct ELISAs and Western blots. In direct ELISAs and Western blots, approximately 10% cross-reactivity with recombinant human MDGA2 is observed.                                                                                                    |
| <b>Source</b>             | Polyclonal Sheep IgG                                                                                                                                                                                                                                                        |
| <b>Purification</b>       | Antigen Affinity-purified                                                                                                                                                                                                                                                   |
| <b>Immunogen</b>          | Mouse myeloma cell line NS0-derived recombinant human MDGA1<br>Gln19-Gly931<br>Accession # Q8NFP4                                                                                                                                                                           |
| <b>Conjugate</b>          | Alexa Fluor 532<br>Excitation Wavelength: 534 nm<br>Emission Wavelength: 553 nm                                                                                                                                                                                             |
| <b>Formulation</b>        | Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide<br><br>*Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions. |

## APPLICATIONS

**Please Note:** Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

**Western Blot** Optimal dilution of this antibody should be experimentally determined.

## PREPARATION AND STORAGE

**Shipping** The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.

**Stability & Storage** Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

## BACKGROUND

MDGA1 (MAM domain containing glycosylphosphatidylinositol anchor 1; also named GPIM) is a 140 kDa member of the Ig superfamily of proteins (1). Human MDGA1 is synthesized as a 955 amino acid (aa) precursor that contains an 18 aa signal sequence, a 914 aa mature chain, and a 23 aa propeptide. The mature chain consists of six Ig-like domains, a fibronectin type III domain (aa 640-739), a MAM domain (aa 751-918), 11 potential sites for N-linked glycosylation, and a GPI anchor (aa 932) (1-3). Alternative splicing in human MDGA1 produces an isoform variant that shows a 48 aa substitution for aa corresponding to positions 926-955. Mature human MDGA1 shares 96%, 94%, and 80% aa sequence identity with mature rat, mouse, and chick MDGA1, respectively. MDGA1 is structurally similar to other IgCAMs, such as the L1 family and axonin 1, which have roles in cell adhesion, migration, and process outgrowth (4). Expression analysis in normal human tissues reveals that MDGA1 is expressed as a 5 kb and 9.5 kb mRNA (2). The smaller transcript is highly expressed in some human cancer cell lines, as well as in different primary tumors (lung, colon, uterus, stomach, and breast) (2). Gene expression is higher in several tumor tissues as compared with corresponding normal tissues, and may be involved in tumor progression (2). The 9.5 kb transcript is found exclusively in the central and peripheral nervous system (1-2, 4). One study shows that MDGA1 is expressed in mice neocortex by layer 2/3 neurons throughout their development (4), and it is involved in controlling the migration and lamination of layer 2/3 neurons (5).

## PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc., and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.