

DESCRIPTION

Species Reactivity	Human
Specificity	Detects endogenous human NEDD9 in Western blots.
Source	Polyclonal Goat IgG
Purification	Antigen Affinity-purified
Immunogen	<i>E. coli</i> -derived recombinant human NEDD9 Lys153-Asn366 Accession # Q14511
Conjugate	Alexa Fluor 405 Excitation Wavelength: 405 nm Emission Wavelength: 421 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide
*Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.	

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

Western Blot Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

NEDD9 (also HEF1 and Cas-L) is a 93 kDa (predicted) member of the CAS family of proteins. It is widely expressed, and serves to integrate integrin β1 and BCR signaling with cell mitotic pathways. Human NEDD9 is 834 amino acids (aa) in length and contains one SH3 domain (aa 3-65), a Tyr-rich SH2-binding region (aa 90-350), one Ser-rich four α-helix bundle protein interaction site (aa 350-650), and an HLH motif that mediates Cas protein dimerization (aa 710-760). NEDD9 usually runs at 105 or 115 kDa in SDS-PAGE due to phosphorylation. It is also cleaved near Leu361 to generate an active N-terminal 55 kDa fragment. One isoform shows a deletion of aa 155-188, while two others show 20 and seven aa substitutions for aa 155-834 and 1-153, respectively. Over aa 153-366, human and mouse NEDD9 share 85% aa identity.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.