

DESCRIPTION


Species Reactivity	Human/Mouse
Specificity	Detects endogenous human/mouse HMGB3 in Western blots.
Source	Polyclonal Goat IgG
Purification	Antigen Affinity-purified
Immunogen	<i>E. coli</i> -derived recombinant human HMGB3 Met1-Val180 Accession # O15347
Formulation	Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose. See Certificate of Analysis for details. *Small pack size (-SP) is supplied either lyophilized or as a 0.2 µm filtered solution in PBS.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

	Recommended Concentration	Sample
Western Blot	1 µg/mL	See Below

DATA



PREPARATION AND STORAGE

Reconstitution	Reconstitute at 0.2 mg/mL in sterile PBS.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below. *Small pack size (-SP) is shipped with polar packs. Upon receipt, store it immediately at -20 to -70 °C
Stability & Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <ul style="list-style-type: none"> 12 months from date of receipt, -20 to -70 °C as supplied. 1 month, 2 to 8 °C under sterile conditions after reconstitution. 6 months, -20 to -70 °C under sterile conditions after reconstitution.

BACKGROUND

High mobility group protein B3 (HMGB3; also HMG-4 and HMG-2a) is a 23 kDa nuclear protein and member of the HMGB family. Human HMGB3 is synthesized as a 200 amino acid (aa) precursor, which is demethylated to produce the 199 aa mature chain. The protein contains two HMG box DNA-binding domains (aa 9-79 and 93-161) and an acidic Asp/Glu-rich region (aa 181-200). Human HMGB3 shares 98% aa sequence identity with mouse and bovine HMGB3. HMGB3 is expressed predominantly in the placenta. It binds preferentially single-stranded DNA and unwinds double stranded DNA.