

DESCRIPTION


Species Reactivity	Mouse
Specificity	Detects mouse TSLP in direct ELISAs and Western blots.
Source	Polyclonal Goat IgG
Purification	Antigen Affinity-purified
Immunogen	<i>S. frugiperda</i> insect ovarian cell line Sf 21-derived recombinant mouse TSLP Tyr20-Glu140 Accession # Q9JIE6
Endotoxin Level	<0.10 EU per 1 µg of the antibody by the LAL method.
Formulation	Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose. See Certificate of Analysis for details. *Small pack size (-SP) is supplied either lyophilized or as a 0.2 µm filtered solution in PBS.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

	Recommended Concentration	Sample
Western Blot	0.1 µg/mL	Recombinant Mouse TSLP (Catalog # 555-TS)
Neutralization		Measured by its ability to neutralize TSLP-induced proliferation in the BaF3 mouse pro-B cell line transfected with mouse IL-7 Rα/CD127 [Park, L.S. <i>et al.</i> (2000) <i>J. Exp. Med.</i> 192 :659]. The Neutralization Dose (ND ₅₀) is typically 0.1-0.4 µg/mL in the presence of 7.5 ng/mL Recombinant Mouse TSLP.

DATA



PREPARATION AND STORAGE

Reconstitution	Reconstitute at 0.2 mg/mL in sterile PBS.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below. *Small pack size (-SP) is shipped with polar packs. Upon receipt, store it immediately at -20 to -70 °C
Stability & Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <ul style="list-style-type: none"> ● 12 months from date of receipt, -20 to -70 °C as supplied. ● 1 month, 2 to 8 °C under sterile conditions after reconstitution. ● 6 months, -20 to -70 °C under sterile conditions after reconstitution.

BACKGROUND

Stromal Lymphopoietin (TSLP) was originally identified from the conditioned medium of a mouse thymic stromal cell line as a protein that promoted the development of B cells. The activity of mouse TSLP overlaps with, but is distinct from, that of mouse IL-7 (1). Mouse TSLP cDNA encodes a 140 amino acid (aa) residue precursor protein with a 19 aa signal sequence. Within the mature region, there are three potential N-glycosylation sites. The Sf 21 cell expressed rmTSLP is likely to be glycosylated at all three sites, as three major glycoforms were visible on SDS-PAGE (Figure 1). Insect cells are known to express relatively simple and homogeneous N-glycans that mainly belong to the high mannose type (2). This recombinant protein was found to be an excellent substrate for N-specific glycosidases such as Endo F3 (Figure 1). The majority of the glycans on rmTSLP can be readily removed by Endo F3. However, a small percentage of the glycans is somewhat resistant to Endo F3 digestion, possibly lacking core fucose, as it is known that core fucosylated N-glycans are strongly preferred substrates for Endo F3 digestion (3).

References:

1. Sims, J.E. *et al.* (2000) *J. Exp. Med.* **192**:671.
2. Staudacher, E. *et al.* (1992) *Eur. J. Biochem.* **207**:987.
3. Tarentino, A.L. and T.H. Jr. Plummer (1994) *Glycobiology* **4**: 771.