

DESCRIPTION

Species Reactivity	Human/Mouse
Specificity	Detects human and mouse EED in direct ELISAs and Western blots.
Source	Polyclonal Sheep IgG
Purification	Antigen Affinity-purified
Immunogen	<i>E. coli</i> -derived recombinant human EED Ser2-Cys150 Accession # O75530
Conjugate	Alexa Fluor 700 Excitation Wavelength: 675-700 nm Emission Wavelength: 723 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide
*Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.	

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

Western Blot Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

EED (Embryonic ectoderm development; also WAIT-1) is a 51 kDa member of the WD repeat ESC family of proteins. It is widely expressed, and constitutes part of both the PRC1 and PRC2 nuclear complexes. PRC1 contains EED, Suz12, RbAp46 and Ezh1 and is found in nonproliferating tissue. PRC2 differs only by the presence of Ezh2 and is found in proliferating tissues. Both complexes methylate histone H3 and act as gene repressors. Human EED is 441 amino acids in length. It contains three phosphorylation sites (Ser34; Thr55; Thr57) followed by seven WD40 domains (aa 91-441). There are two alternate start sites, one at Met15, and another that is 94 aa upstream of the standard site, begins with a Val and generates a 69 kDa isoform. A third isoform variant shows an insertion of 25 aa after Ser323, while a fourth variant shows a premature truncation after Lys400. Human and mouse EED are identical in amino acid sequence.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.