

DESCRIPTION

Species Reactivity	Human
Specificity	Detects human AIRE in direct ELISAs and Western blots.
Source	Polyclonal Goat IgG
Purification	Antigen Affinity-purified
Immunogen	<i>E. coli</i> -derived recombinant human AIRE Ser476-Ser545 Accession # O43918
Conjugate	Alexa Fluor 488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide
*Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.	

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

Western Blot Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

AIRE (Autoimmune Regulator; also APECED) is a 58-60 kDa member of the PHD Zn-finger family of transcription factors. It appears to have restricted expression, being identified in thymic medullary epithelium, dendritic cells and neutrophils. AIRE regulates the expression of self-antigens in thymus. In particular, genes coding for these antigens are generally inactive, and show minimal H3K4me3 and AcH3 reactivity in their promoters. Human AIRE is 545 amino acids (aa) in length. It contains multiple domains, including an HSR domain that mediates homodimerization (aa 1-105), an NLS (aa 131-133), a SAND domain (aa 181-280), and two Zn-finger PHD domains (aa 296-343 and 434-475) that bind to nonmethylated histone H3. There are isoforms that show a Trp substitution for aa 437-467, a 95 aa substitution for aa 1-292, and a 75 aa substitution for aa 377-545. Over aa 476-545, human AIRE shares 71% aa identity with mouse AIRE.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.