

DESCRIPTION

Species Reactivity	Human/Mouse/Rat
Specificity	Detects human, mouse, and rat SNAP25 in direct ELISAs and Western blots. In direct ELISAs, approximately 5% cross-reactivity with recombinant human SNAP23 is observed.
Source	Polyclonal Sheep IgG
Purification	Antigen Affinity-purified
Immunogen	<i>E. coli</i> -derived recombinant human SNAP25 Asp140-Gly206 Accession # P60880
Conjugate	Alexa Fluor Plus 488 Excitation Wavelength: 493 nm Emission Wavelength: 518 nm
Formulation	Supplied 0.2 mg/mL in a saline solution containing BSA and Sodium Azide. *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. General Protocols are available in the Technical Information section on our website.

Western Blot	Optimal dilution of this antibody should be experimentally determined.
Immunocytochemistry	Optimal dilution of this antibody should be experimentally determined.

DATA

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

SNAP25 (Synaptosomal associated protein of 25 kDa; also SUP) is a cytosolic 24-29 kDa member of the SNAP25 family of proteins. It is expressed in neurons, β-cells, and adrenal chromaffin cells. SNAP25 is involved in vesicle exocytosis, where two presynaptic SNAP25, and one syntaxin-1 molecule associates with vesicle membrane VAMP2 to form a α-helix via a zipper-like mechanism. This approximates two membranes, which subsequently fuse to create a pore. Human SNAP25b is 206 amino acids (aa) in length. It contains two t-SNARE coiled-coil homology domains (aa 19-81 and 140-202) and a proteolytic cleavage site at Arg180:Ile181. Palmitoylation occurs between aa 85-92; phosphorylation occurs on Thr138. There are three isoform variants. One shows nine aa changes over aa 58-89 (isoform a), a second shows a 30 aa substitution for aa 1-94, while a third contains an alternative start site at Met64. Full-length human and mouse SNAP25 are identical in amino acid sequence.

PRODUCT SPECIFIC NOTICES

This product is provided under an intellectual property license from Life Technologies Corporation. The transfer of this product is conditioned on the buyer using the purchased product solely in research conducted by the buyer, excluding contract research or any fee for service research, and the buyer must not (1) use this product or its components for (a) diagnostic, therapeutic or prophylactic purposes; (b) testing, analysis or screening services, or information in return for compensation on a per-test basis; or (c) manufacturing or quality assurance or quality control, and/or (2) sell or transfer this product or its components for resale, whether or not resold for use in research. For information on purchasing a license to this product for purposes other than as described above, contact Life Technologies Corporation, 5781 Van Allen Way, Carlsbad, CA 92008 USA or outlicensing@thermofisher.com.