

DESCRIPTION

Species Reactivity	Human/Mouse
Specificity	Detects human Proprotein Convertase 2/PCSK2 in direct ELISAs and, human and mouse Proprotein Convertase 2/PCSK2 in Western blots.
Source	Polyclonal Sheep IgG
Purification	Antigen Affinity-purified
Immunogen	Chinese hamster ovary cell line CHO-derived recombinant human Proprotein Convertase 2/PCSK2 Gly110-Asn638 Accession # P16519
Conjugate	Alexa Fluor 532 Excitation Wavelength: 534 nm Emission Wavelength: 553 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide
*Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.	

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the Technical Information section on our website.

Western Blot	Optimal dilution of this antibody should be experimentally determined.
Immunocytochemistry	Optimal dilution of this antibody should be experimentally determined.
Immunoprecipitation	Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

Proprotein Convertase 2 (PCSK2) is a subtilisin-like serine peptidase that processes proteins into biologically active products in the secretory pathway. Like other furin/kexin family members, PCSK2 cleaves its substrates primarily after sites of paired basic amino acid residues (1). Autoactivation of PCSK2 occurs in a post-Golgi compartment of the secretory system (2). Interaction with the secretory protein 7B2 is required for proper activation of PCSK2 (3). The N-terminal domain of 7B2 stabilizes active PCSK2, and a C-terminal fragment of 7B2 can inhibit PCSK2 activity (4). rhPCSK2 was coexpressed with 7B2 to facilitate activation of the enzyme. PCSK2 is a major proteolytic processing enzyme of the regulated secretory pathway of the neuroendocrine system, where it generates a number of hormones and neuropeptides (5). Products of PCSK2 processing include enkephalins, insulin, somatostatin, dynorphin, and LHRH. Within the region used as an immunogen, human PCSK2 shares 97% aa identity with mouse PCSK2.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.