

DESCRIPTION

Species Reactivity	Human
Specificity	Detects human HGK in Western blots.
Source	Polyclonal Goat IgG
Purification	Antigen Affinity-purified
Immunogen	<i>E. coli</i> -derived recombinant human HGK Ser686-Val820 Accession # O95819-4
Conjugate	Alexa Fluor Plus 594 Excitation Wavelength: 590 nm Emission Wavelength: 618 nm
Formulation	Supplied 0.2 mg/mL in a saline solution containing BSA and Sodium Azide. *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. General Protocols are available in the Technical Information section on our website.

Western Blot Optimal dilution of this antibody should be experimentally determined.

DATA

PREPARATION AND STORAGE

Shipping The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.

Stability & Storage Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

HGK (Hepatocyte kinase-like/Germinal center like-kinase Kinase; also MAP4K4) is a 140 kDa member of the MSN subfamily, MAP4K family of molecules. It is ubiquitously expressed, and intimately associated with TNF α signaling. HGK is activated by TNF-induced c-Jun activation, and blocks both PPAR γ expression and glucose-stimulated insulin secretion. Human HGK is 1239 amino acids (aa) in length (Accession # O95819) and contains one protein kinase domain (aa 25-289) and a CHN region (aa 926-1213). There are multiple splice variants. The antibody product offered by R&D Systems uses a 130 kDa, 1166 aa isoform as immunogen. It shows a deletion of aa 495-525 and 570-623, plus an eight aa insertion after His1112, and a GEV insertion after Ala740.

PRODUCT SPECIFIC NOTICES

This product is provided under an intellectual property license from Life Technologies Corporation. The transfer of this product is conditioned on the buyer using the purchased product solely in research conducted by the buyer, excluding contract research or any fee for service research, and the buyer must not (1) use this product or its components for (a) diagnostic, therapeutic or prophylactic purposes; (b) testing, analysis or screening services, or information in return for compensation on a per-test basis; or (c) manufacturing or quality assurance or quality control, and/or (2) sell or transfer this product or its components for resale, whether or not resold for use in research. For information on purchasing a license to this product for purposes other than as described above, contact Life Technologies Corporation, 5781 Van Allen Way, Carlsbad, CA 92008 USA or outlicensing@thermofisher.com.