

DESCRIPTION

Species Reactivity	Human/Mouse
Specificity	Detects human and mouse SIN3A in Western blots.
Source	Polyclonal Goat IgG
Purification	Antigen Affinity-purified
Immunogen	<i>E. coli</i> -derived recombinant human SIN3A Tyr1187-Pro1273 Accession # Q96ST3
Conjugate	Alexa Fluor 750 Excitation Wavelength: 749 nm Emission Wavelength: 775 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide
*Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.	

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

Western Blot Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

SIN3A (Switch insensitive-3 A) is a nuclear, 140-150 kDa member of the SIN3 family of co-repressors. SIN3A plays a key role in gene repression. Normally, genes are not expressed when chromatin is "closed", or condensed, precluding transcriptional activators from assembling on gene promoters. This state is achieved through histone deacetylation and DNA methylation, and SIN3A serves as a scaffold for multiple proteins, including deacetylases (HDAC1 and 2). Human SIN3A is 1273 amino acids (aa) in length. It contains three PAH (paired amphipathic helix) repeats (aa 141-187, 322-381 and 477-523) and multiple phosphoserines and acetylated lysines. There are three potential splice variants. The first shows a 33 aa substitution for aa 124-1273, while the second and third show eight aa substitutions for aa 337-352, coupled to either a deletion of aa 1009-1066, or a two aa substitution for aa 1008-1065. Over aa 1187-1273, human and mouse SIN3A show 100% identity in aa sequence.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.