

DESCRIPTION

Species Reactivity	Mouse
Specificity	Detects mouse MEPE/OF45 in direct ELISAs and Western blots. In direct ELISAs, approximately 10% cross-reactivity with recombinant human MEPE/OF45 is observed.
Source	Polyclonal Sheep IgG
Purification	Antigen Affinity-purified
Immunogen	Mouse myeloma cell line NS0-derived recombinant mouse MEPE/OF45 Ala17-Asp433 Accession # AAK70342
Conjugate	Alexa Fluor 488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

Western Blot Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.

Stability & Storage Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

MEPE (matrix extracellular phosphoglycoprotein; also OF45) is a 56-58 kDa member of the SIBLING protein family. It is expressed by osteocytes, osteoblasts, odontoblasts and chondrocytes. MEPE is found both intra- and extracellularly. Within the cell, it binds to CHK1. This protects CHK1 from degradation, and the cell from DNA damage. Extracellularly, it binds to Phex, which protects it from cathepsin B proteolysis. Undegraded, it promotes bone vascularization and blocks bone resorption. Following cleavage after Arg410 (in mouse), a 2 kDa ASARM fragment is generated that, when phosphorylated, blocks mineralization. Mature mouse MEPE is 417 amino acids (aa) in length. It contains an RGD motif and ASARM domain (aa 416-433). In addition to its a 57 kDa form, 70 kDa, 42-45 kDa and 30-35 kDa forms of mouse MEPE have been reported. Whether they are cleavage forms or products of glycosylation is unclear. Over aa 147-433, mouse MEPE shares 77% and 57% aa identity with rat and human MEPE, respectively.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.