

DESCRIPTION

Species Reactivity	Human
Specificity	Detects human Leucyl-cystinyl Aminopeptidase/LNPEP in direct ELISAs and Western blots. In direct ELISAs, less than 2% cross-reactivity with recombinant human (rh) ENPEP, rhLRAP, and rhARTS1 is observed.
Source	Polyclonal Sheep IgG
Purification	Antigen Affinity-purified
Immunogen	Chinese hamster ovary cell line CHO-derived recombinant human Leucyl-cystinyl Aminopeptidase/LNPEP Arg133-Gly552 Accession # Q9UIQ6
Conjugate	Alexa Fluor 488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

Western Blot	Optimal dilution of this antibody should be experimentally determined.
Immunoprecipitation	Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

LNPEP (also known as P-LAP, CAP, IRAP, and Oxytocinase) is a 165-175 kDa member of the oxytocinase subfamily, M1 family of zinc metallopeptidases. It is widely expressed after cell stimulation and is found on endothelial cells, trophoblasts, skeletal muscle cells, neurons, adipocytes, and renal distal tubule epithelium. LNPEP cleaves vasopressin and oxytocin and serves as a receptor for AngIV. Human LNPEP is a 1025 amino acid (aa) type II transmembrane glycoprotein (SwissProt #:Q9UIQ6). It contains a 110 aa cytoplasmic region (aa 1-110) plus an 894 aa extracellular domain (aa 132-1025) that contains a peptidase M1 region (aa 168-552). There are multiple potential splice variants. There are two alternative start sites at Met15 and Met20, and the cytoplasmic tail undergoes phosphorylation at Tyr70 and Ser80. LNPEP is reportedly forms homodimers. Proteolytic cleavage between Phe154Ala155 generates a 150-155 kDa soluble form. Over aa 133-552, human LNPEP shares 89% aa sequence identity with mouse LNPEP.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.