

DESCRIPTION

Species Reactivity	Human
Specificity	Detects human KDM5B in Western blots.
Source	Polyclonal Sheep IgG
Purification	Antigen Affinity-purified
Immunogen	<i>E. coli</i> -derived recombinant human KDM5B Leu819-Arg918 Accession # Q9UGL1
Conjugate	Alexa Fluor 647 Excitation Wavelength: 650 nm Emission Wavelength: 668 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide
*Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.	

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. General Protocols are available in the Technical Information section on our website.

Western Blot	Optimal dilution of this antibody should be experimentally determined.
Immunocytochemistry	Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

KDM5B (also JARID1B, RBP2-H1 and PLU-1) is a 170-200 kDa member of the JARID1 histone demethylase family of proteins. It is a transcriptional repressor that is highly expressed in Sertoli cells and mammary epithelium during pregnancy, and interacts with BF-1 and PAX9. KDM5B is reported to demethylate tri- and di-methylated Lys4 on histone H3. It has no activity on Lys9 or 27. Human KDM5B is 1544 amino acids (aa) in length. It contains a jumonji N domain (aa 32-73), an ARID DNA binding region (aa 97-187), a PHD domain (aa 309-359), a jumonji C domain (aa 486-602), one C5HC2 motif (aa 692-745), and two consecutive PHD domains (aa 1176-1538). There are three utilized phosphorylation sites at Ser1383, 1384 and 1405. There are three potential splice variants. Two show a premature truncation after Pro1433 coupled to an alternative start site at either Met396 or Met159. A third contains a 36 aa insertion after Glu237. Over aa 819-918, human KDM5B shares 96% aa identity with mouse KDM5B.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.