

DESCRIPTION

Species Reactivity	Human
Specificity	Detects human ASCL2/Mash2 in Western blots.
Source	Polyclonal Sheep IgG
Purification	Antigen Affinity-purified
Immunogen	<i>E. coli</i> -derived recombinant human ASCL2/Mash2 Met1-Ala49 Accession # Q99929
Conjugate	Alexa Fluor 750 Excitation Wavelength: 749 nm Emission Wavelength: 775 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide
*Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.	

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

Western Blot	Optimal dilution of this antibody should be experimentally determined.
Immunohistochemistry	Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

ASCL2 (Achaete-scute like 2; also HASH2, bHLHa45 and Mash2) is a class II member of the HLH family of transcription factors. Its predicted MW is 20 kDa. ASCL2 shows restricted expression, being limited to intestinal Lgr5⁺ stem cells and first trimester placental cytotrophoblasts. In the intestine, ASCL2 is under the control of Wnts and serves to maintain an epithelium stem cell pool. In the placenta, ASCL2 acts to maintain the pool of cytotrophoblasts at the expense of syncytiotrophoblasts, thus promoting placental growth. Human ASCL2 is 193 amino acids (aa) in length. It contains one poly-Arg motif (aa 36-39), a DNA binding sequence (aa 53-63) and an HLH domain (aa 64-103). ASCL2 forms heterodimers with daughterless homologs (E22, E2A and HEB). Over aa 1-49, human ASCL2 shares 59% aa identity with mouse Mash2.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.