

DESCRIPTION

Species Reactivity	Human
Specificity	Detects human SUMO1-Specific Peptidase 1/SEN1 in direct ELISAs.
Source	Polyclonal Sheep IgG
Purification	Antigen Affinity-purified
Immunogen	<i>E. coli</i> -derived recombinant human SUMO1-Specific Peptidase 1/SEN1 Glu419-Leu644 Accession # Q9P0U3
Conjugate	Alexa Fluor 405 Excitation Wavelength: 405 nm Emission Wavelength: 421 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide
*Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.	

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

Immunocytochemistry	Optimal dilution of this antibody should be experimentally determined.
Immunohistochemistry	Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

SEN1 (Sentrin/SUMO-specific protease) is a member of the SENP family of proteases. Sentrin/SUMO-specific proteases (SENPs) are a group of cysteine-type peptidases that catalyze two essential functions in the SUMO pathways: processing of full-length SUMOs to their mature forms and deconjugation of SUMOs from SUMOylated proteins. The seven mammalian SENPs share a conserved C-terminal catalytic domain while the N-terminal domains have no significant similarity. Human SENP-1 has broad specificity for the three mammalian SUMOs. It is found in the cytoplasm and nucleus depending on cell type, and is expressed in testis, thymus, pancreas, spleen, liver, ovary and small intestine. It is thought that localization of SENP-1 is vital for the regulation for SUMOylation status of target proteins. The small ubiquitin-like modifier (SUMO) is a member of ubiquitin-like protein family. SUMO modification of the target proteins is a reversible process that regulates many cellular processes including transcription regulation, nuclear localization, centromere segregation and signal transduction. The recombinant human SENP-1 (immunogen) includes the catalytic domain, which has been shown to be sufficient for SENP-1 activity and substrate specificity.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.