

DESCRIPTION

Species Reactivity	Human/Mouse/Rat
Specificity	Detects recombinant mouse and rat GABA _B R1 in direct ELISAs and Western blots. Detects human and rat GABA _B R1 in Western blots. In direct ELISAs, less than 1% cross-reactivity with recombinant rat GABA _B R2 is observed.
Source	Polyclonal Sheep IgG
Purification	Antigen Affinity-purified
Immunogen	Chinese hamster ovary cell line CHO-derived recombinant rat GABA _B R1 Gly17-Leu586 Accession # Q9Z0U4
Conjugate	Alexa Fluor 700 Excitation Wavelength: 675-700 nm Emission Wavelength: 723 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the Technical Information section on our website.

Western Blot	Optimal dilution of this antibody should be experimentally determined.
Immunohistochemistry	Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

GABA_B R1 (GABA-B receptor subunit 1; also GABA-BR1, GABBR1 and GB1) is a multispan member of the GABA-B receptor subfamily, GPCR-3 family of proteins. It forms an obligatory heterodimer with GABA-BR2, creating a G-protein metabotropic GABA receptor that inhibits adenylyl cyclase activity and activates K⁺ channels. Presynaptically, this blocks neurotransmitter release; postsynaptically, it lowers neuron excitability. Rat GABA_B R1 is 991 amino acids (aa) in length. It is a 7-transmembrane glycoprotein that contains a 16 aa signal sequence, an extended N-terminal extracellular region (aa 17-590) that contains two SUSHI domains (aa 29-158), and a long C-terminal cytoplasmic domain (aa 885-991). There are several splice variants with predicted molecular weights ranging from 90 to 111 kDa and multiple glycosylation sites. The 991 aa isoform described above is called GABA_B R1e (R1e). There is also a 960 aa, 130 kDa isoform that shows a deletion of aa 771-801. This variant (R1a) is associated with postsynaptic membranes. A third isoform (R1b) is 844 aa in length and 100 kDa in size, and possesses both a deletion of aa 771-801, and a 47 aa substitution for aa 1-163. This variant is presynaptic in location. Two other isoforms are variants of GABA_B R1b. Each show the same N-terminal substitution, with a fourth isoform (R1c) retaining aa 771-801, and a fifth isoform (R1d) deleting aa 771-801, coupled to a 25 aa substitution for aa 935-991. Over aa 17-586, rat GABA_B R1e/a shares 99% aa identity with both mouse and human GABA_B R1.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc., and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.