

DESCRIPTION

Species Reactivity	Mouse
Specificity	Detects mouse NELL1 in direct ELISAs and Western blots. In direct ELISAs, approximately 25% cross-reactivity with recombinant human NELL1 is observed and approximately 10% cross-reactivity with recombinant mouse NELL2 is observed.
Source	Polyclonal Sheep IgG
Purification	Antigen Affinity-purified
Immunogen	Chinese hamster ovary cell line CHO-derived recombinant mouse NELL1 Arg17-Asn810 Accession # Q2VWQ2
Conjugate	Alexa Fluor 700 Excitation Wavelength: 675-700 nm Emission Wavelength: 723 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

Western Blot	Optimal dilution of this antibody should be experimentally determined.
Immunohistochemistry	Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

NELL1 (neural EGF-like like protein 1) is an 89 kDa (predicted) member of the EGF-like domain containing family, Laminin G/N-TSP1/Pentraxin gene superfamily of molecules. When secreted, NELL1 exists as a phosphoglycoprotein that can add as much as 50 kDa to the calculated MW. NELL1 has restricted expression, being limited to pre-B cells and osteoblasts, where it apparently promotes osteoblast maturation and bone formation. In tumors, it is found in neuroblastoma-derived cells. NELL1 is both secreted and retained intracellularly where it is phosphorylated by PKC. The mouse NELL1 precursor is 810 amino acids (aa) in length. It contains a 16 aa signal sequence plus a 794 aa mature region. The mature region possesses an N-terminal TSP domain (aa 81-230), two VWFC domains (aa 271-390), six consecutive EGF-like domains (aa 391-631), and three additional C-terminal VWFC domains (aa 632-807). Secreted NELL1 forms a 400-420 kDa noncovalent homotrimer. Over aa #17-810, mouse NELL1 shares 98% and 93% aa identity with rat and human NELL1, respectively.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc., and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.