

DESCRIPTION

Species Reactivity	Human/Mouse
Specificity	Detects human and mouse MyBPC3 in Western blots and in direct ELISAs.
Source	Polyclonal Sheep IgG
Purification	Antigen Affinity-purified
Immunogen	<i>E. coli</i> -derived recombinant mouse MyBPC3 Arg998-Lys1100 Accession # O70468
Conjugate	Alexa Fluor 532 Excitation Wavelength: 534 nm Emission Wavelength: 553 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. General Protocols are available in the Technical Information section on our website.

Western Blot Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

MyBPC3 (Myosin-binding protein C-cardiac type) is a 140-150 kDa member of the MyBP family, Ig superfamily of proteins. It is expressed in cardiac muscle, and contributes both to myosin filament structure by interacting with light meromyosin, and the regulation of contraction by binding to myosin subfragment-2, which results in a reduction of actomyosin ATPase activity. Mouse MyBPC3 is 1270 amino acids (aa) in length. It contains five consecutive C2-type Ig-like domains (aa 151-767), three FN type III repeats (aa 768-958), and two additional C-terminal Ig-like domains (aa 967-1270). There are at least three utilized phosphorylation sites and one Pro-rich region (aa 100-150). In human, multiple mutations generate variable-length premature truncated forms of MyBPC3. There are two potential isoform variants. One basically shows a six aa substitution for aa 339-344, while a second contains the same six aa substitution coupled to an eight aa extension at the N-terminus. Over aa 998-1100, mouse MyBPC3 shares 89% and 95% aa identity with human and rat MyBPC3, respectively.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.