

DESCRIPTION

Species Reactivity	Mouse
Specificity	Detects mouse SIGNR3/CD209d in direct ELISAs and Western blots. In direct ELISAs, less than 1% cross-reactivity with recombinant human (rh) DC-SIGN R, rhDC-SIGN, recombinant mouse (rm) SIGN R1, rmSIGN R4, and rmSIGN R7 is observed.
Source	Polyclonal Sheep IgG
Purification	Antigen Affinity-purified
Immunogen	Mouse myeloma cell line NS0-derived recombinant mouse SIGNR3/CD209d Glu82-Lys237 Accession # Q91ZW8
Conjugate	Alexa Fluor 488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

Western Blot Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.

Stability & Storage Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

SIGNR3 (Specific Intercellular adhesion molecule-3-Grabbing Non-integrin Related gene 3; also CD209d and DC-SIGNR3) is a 30-34 kDa member of the C-type lectin family of proteins. It is found on macrophages, DC and monocytes, and appears to bind both high mannose and fucose moieties on the surface of extracellular pathogens. It is particularly important in clearing Mycobacterium tuberculosis. Upon ligation, SIGNR3 is internalized and presumably enters the endosomal system. Notably, and unlike in human, there are multiple SIGNs in mouse, and SIGNR3 most closely resembles the biology of human DC-SIGN. Mouse SIGNR3 is a 237 amino acid (aa) type II transmembrane glycoprotein. It contains a 54 aa N-terminal cytoplasmic region plus a 162 aa C-terminal extracellular domain (aa 76-237) that possesses a C-type lectin domain (aa 106-228). It is suggested that mouse SIGNR3 may oligomerize. There are two potential splice variants. One possesses an alternative start site 23 aa upstream of the standard site, while a second shows a deletion of aa 45-74. Over aa 82-237, mouse SIGNR3 shares 89% and 59% aa sequence identity with rat and human SIGNR3, respectively.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.