

DESCRIPTION

Species Reactivity	Human
Specificity	Detects human ST6 Gal Sialyltransferase 2/ST6GAL2 in direct ELISAs and Western blots. In direct ELISAs, approximately 10% cross-reactivity with recombinant human ST6GAL1 is observed.
Source	Polyclonal Sheep IgG
Purification	Antigen Affinity-purified
Immunogen	Chinese hamster ovary cell line CHO-derived recombinant human ST6 Gal Sialyltransferase 2/ST6GAL2 Glu243-Ser529 Accession # Q96JF0
Conjugate	Alexa Fluor 594 Excitation Wavelength: 590 nm Emission Wavelength: 617 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

Western Blot	Optimal dilution of this antibody should be experimentally determined.
Immunohistochemistry	Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

ST6GAL2 (Sialyltransferase 2; also β -galactoside α -2,6 Sialyltransferase 2) is a 60-62 kDa member of the glycosyltransferase 29 family of enzymes. It has restricted expression, likely being found in neurons and bronchial epithelium. ST6GAL2 transfers sialic acid from CMP-Neu5Ac to a galactose that is linked to N-Acetylglucosamine through a β 1-4 linkage (Gal β 1-4GlcNAc). This transfer will occur when the target motif is found either on select glycoproteins or free oligosaccharides. Notably, the kinetics of this enzyme suggest additional biological functions other than sialic acid transfer. Human ST6GAL2 is a 529 amino acid (aa) type II transmembrane Golgi protein. It contains an 11 aa N-terminal cytoplasmic region and a 497 aa luminal domain (aa 33-529). The luminal domain possesses an L motif that binds CMP-Sialic acid (aa 296-343), an S motif that binds CMP-Sialic acid and galactose (aa 444-458) and a VS motif that is involved with the enzyme's catalytic activity (aa 483-488). There is some evidence that ST6GAL2 may form homodimers. There is one splice variant that shows a 26 aa substitution for aa 441-529. Over aa 243-529, human ST6GAL2 shares 88% aa sequence identity with mouse ST6GAL2.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.