

DESCRIPTION

Species Reactivity	Human
Specificity	Detects human USE1/UBE2Z in direct ELISAs and Western blots.
Source	Polyclonal Sheep IgG
Purification	Antigen Affinity-purified
Immunogen	<i>E. coli</i> -derived recombinant human Use1/UBE2Z Met1-Val354 Accession # Q9H832
Conjugate	Alexa Fluor 700 Excitation Wavelength: 675-700 nm Emission Wavelength: 723 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide
*Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.	

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

Western Blot Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

Uba6-specific E2 Conjugating Enzyme 1 (Use1), also known as Ubiquitin-conjugating Enzyme E2Z (UBE2Z), is a 38 to 43 kDa member of the Ubiquitin-conjugating (E2) enzyme family that is highly expressed in human placenta, pancreas, spleen, and testis. Use1 has an E2 catalytic core domain that contains an active site cysteine residue that is required for the formation of a thioester bond with Ubiquitin. It functions specifically with the UBE1L2/UBA6 Ubiquitin-activating (E1) enzyme, rather than the UBE1 E1 enzyme, for Ubiquitin conjugation. Expression is enriched in the cytoplasm, where it mediates the ubiquitination of RGS4 and RGS5 via the N-end rule proteolytic pathway. Use1 also functions with UBE1L2/UBA6 to mediate the conjugation of HLA F-associated Transcript 10 (FAT10), a Ubiquitin-like protein, and is auto-FAT10ylated in response to TNF-alpha and IFN-gamma. Full-length human Use1 is 354 amino acids (aa) and shares 100% sequence identity with mouse and rat Use1. A potential smaller isoform lacks aa 1-108.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc., and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.