

DESCRIPTION

Species Reactivity	Human
Specificity	Detects human APRIL/TNFSF13 in direct ELISAs and Western blots. In direct ELISAs, less than 10% cross-reactivity with recombinant human (rh) Fas Ligand, rhTNF-α, rhTRAIL, and rhTRANCE is observed.
Source	Polyclonal Goat IgG
Purification	Antigen Affinity-purified
Immunogen	<i>E. coli</i> -derived recombinant human APRIL/TNFSF13
Conjugate	Alexa Fluor 700 Excitation Wavelength: 675-700 nm Emission Wavelength: 723 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

CyTOF-ready	Optimal dilution of this antibody should be experimentally determined.
Western Blot	Optimal dilution of this antibody should be experimentally determined.
ELISA	Optimal dilution of this antibody should be experimentally determined.
Flow Cytometry	Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

APRIL (a proliferation-inducing ligand), also known as TNFSF13, TALL2, and TRDL1, is a member of the TNF ligand superfamily (1). APRIL is synthesized as a 32 kDa type II transmembrane protein which is cleaved by furin in the Golgi to release a 17 kDa soluble molecule (2, 3). Secreted APRIL consists almost entirely of a single TNF homology domain (2, 3). Little or no transmembrane APRIL is expressed on the cell surface (3). Alternate splicing generates isoforms with short deletions at the N- or C-terminus (4). Human APRIL shares 85% aa sequence identity with mouse and rat APRIL. Among TNF superfamily ligands, BAFF shows the greatest sequence homology with APRIL, and the two proteins exhibit overlapping biological activities. APRIL promotes cellular proliferation and protects from apoptosis in normal and transformed cells (3, 5-7). It is present in elevated amounts in a wide variety of cancers primarily due to expression by tumor-infiltrating neutrophils (4, 5, 7-9). Both APRIL and BAFF bind and signal through the TNF superfamily receptors TACI and BCMA, and BAFF additionally functions through BAFF R (6, 10, 11). A stretch of basic amino acids at the N-terminus of APRIL is required for its interaction with heparan sulfate proteoglycans (HSPGs) (12, 13). Binding to HSPGs is independent of APRIL's binding to TACI and BCMA (12, 13). Interaction with HSPGs serves to concentrate APRIL on the surface of cells, thereby favoring TACI- or BCMA-mediated effects (8, 9, 13). APRIL can form bioactive heterotrimers with BAFF, and these circulate in the serum of patients with rheumatic immune disorders (14). A bioactive protein known as TWE-PRIL consists of the intracellular domain, transmembrane segment, and stalk region of TWEAK fused to the TNF homology domain of APRIL (15). TWE-PRIL is expressed in monocytes and activated T cells, and in contrast to APRIL, it is presented on the cell surface (15).

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.