

#### DESCRIPTION




<b>Species Reactivity</b>	Human
<b>Specificity</b>	Detects human SEC13 in direct ELISA and Western blots.
<b>Source</b>	Polyclonal Rabbit IgG
<b>Purification</b>	Antigen Affinity-purified
<b>Immunogen</b>	<i>E. coli</i> -derived recombinant human SEC13 Gly123-Gln322 Accession # P55735
<b>Formulation</b>	Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose. See Certificate of Analysis for details. *Small pack size (-SP) is supplied as a 0.2 µm filtered solution in PBS.

#### APPLICATIONS

**Please Note:** Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

	<b>Recommended Concentration</b>	<b>Sample</b>
<b>Western Blot</b>	0.05 µg/mL	See Below
<b>Immunocytochemistry</b>	1-15 µg/mL	See Below
<b>Simple Western</b>	0.5 µg/mL	See Below

#### DATA



#### PREPARATION AND STORAGE

<b>Reconstitution</b>	Reconstitute at 0.2 mg/mL in sterile PBS.
<b>Shipping</b>	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below. *Small pack size (-SP) is shipped with polar packs. Upon receipt, store it immediately at -20 to -70 °C
<b>Stability &amp; Storage</b>	<b>Use a manual defrost freezer and avoid repeated freeze-thaw cycles.</b> <ul style="list-style-type: none"> <li>• 12 months from date of receipt, -20 to -70 °C as supplied.</li> <li>• 1 month, 2 to 8 °C under sterile conditions after reconstitution.</li> <li>• 6 months, -20 to -70 °C under sterile conditions after reconstitution.</li> </ul>

#### BACKGROUND

Sec13 is a ubiquitously expressed protein which participates in the formation of vesicles in the COPII complex with Sec23p-Sec24p, Sar1p and Sec31. The COPII vesicle budding complex was first described in the yeast system and was shown to provide coat proteins in the early secretory pathway. Sec13 is recruited to endoplasmic reticulum (ER) membranes where it forms a three dimensional cuboctahedron cage lattice structure by association with Sec31. This vesicular structure is necessary for ER-Golgi transport. Presenilin-1, which is implicated in early-onset Alzheimer's disease has been shown to interact with Sec13 vesicles. Additionally, Sec13 has been shown to be a component of the Y-shaped Nup107-160 subcomplex of the nuclear pore complex (NPC), and the GATOR2 sub-complex of the GATOR complex.