

DESCRIPTION

Source	<i>E. coli</i> -derived human TNF-alpha protein Val77-Leu233, with and without an N-terminal Met Accession # P01375 Produced using non-animal reagents in an animal-free laboratory.
N-terminal Sequence Analysis	Met & Val77
Predicted Molecular Mass	17 kDa

SPECIFICATIONS

SDS-PAGE	17 kDa, reducing conditions
Activity	Measured in a cytotoxicity assay using L-929 mouse fibroblast cells in the presence of the metabolic inhibitor actinomycin D. Matthews, N. and M.L. Neale (1987) in <i>Lymphokines and Interferons, A Practical Approach</i> . Clemens, M.J. <i>et al.</i> (eds); IRL Press. 221. The ED ₅₀ for this effect is 25-100 pg/mL. The specific activity of Recombinant Human TNF-α is approximately 7.6 x 10 ⁴ IU/μg, which is calibrated against human TNF-α WHO International Standard (NIBSC code: 88/786).
Endotoxin Level	<0.10 EU per 1 μg of the protein by the LAL method.
Purity	>97%, by SDS-PAGE visualized with Silver Staining and quantitative densitometry by Coomassie® Blue Staining.
Formulation	Lyophilized from a 0.2 μm filtered solution in PBS. See Certificate of Analysis for details.

PREPARATION AND STORAGE

Reconstitution	Reconstitute at 0.2 mg/mL in sterile PBS.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <ul style="list-style-type: none"> • 12 months from date of receipt, -20 to -70 °C as supplied. • 1 month, 2 to 8 °C under sterile conditions after reconstitution. • 3 months, -20 to -70 °C under sterile conditions after reconstitution.

BACKGROUND

Tumor necrosis factor alpha (TNF-α), also known as cachectin and TNFSF2, is the prototypic ligand of the TNF superfamily. It is a pleiotropic molecule that plays a central role in inflammation, immune system development, apoptosis, and lipid metabolism (1, 2). Human TNF-α consists of a 35 amino acid (aa) cytoplasmic domain, a 21 aa transmembrane segment, and a 177 aa extracellular domain (ECD) (3). Within the ECD, human TNF-α shares 97% aa sequence identity with rhesus and 71%-92% with bovine, canine, cotton rat, equine, feline, mouse, porcine, and rat TNF-α. TNF-α is produced by a wide variety of immune, epithelial, endothelial, and tumor cells (1, 2). TNF-α is assembled intracellularly to form a noncovalently linked homotrimer which is expressed on the cell surface (4). Cell surface TNF-α can induce the lysis of neighboring tumor cells and virus infected cells, and it can generate its own downstream cell signaling following ligation by soluble TNFR I (2, 5). Shedding of membrane bound TNF-α by TACE/ADAM17 releases the bioactive cytokine, a 55 kDa soluble trimer of the TNF-α extracellular domain (6-8). TNF-α binds the ubiquitous 55-60 kDa TNF RI (9, 10) and the hematopoietic cell-restricted 80 kDa TNF RII (11, 12), both of which are also expressed as homotrimers (1, 2, 13). Both type I and type II receptors bind TNF-α with comparable affinity (14), although only TNF RI contains a cytoplasmic death domain which triggers the activation of apoptosis. Soluble forms of both types of receptors are released and can neutralize the biological activity of TNF-α (15).

References:

1. Zelova, H. and J. Hosek (2013) *Inflamm. Res.* **62**:641.
2. Juhasz, K. *et al.* (2013) *Expert Rev. Clin. Immunol.* **9**:335.
3. Pennica, D. *et al.* (1984) *Nature* **312**:724.
4. Tang, P. *et al.* (1996) *Biochemistry* **35**:8216.
5. Perez, C. *et al.* (1990) *Cell* **63**:251.
6. Black, R.A. *et al.* (1997) *Nature* **385**:729.
7. Moss, M.L. *et al.* (1997) *Nature* **385**:733.
8. Gearing, A.J.H. *et al.* (1994) *Nature* **370**:555.
9. Schall, T.J. *et al.* (1990) *Cell* **61**:361.
10. Loetscher, H. *et al.* (1990) *Cell* **61**:351.
11. Dembic, Z. *et al.* (1990) *Cytokine* **2**:231.
12. Smith, C.A. *et al.* (1990) *Science* **248**:1019.
13. Loetscher, H. *et al.* (1991) *J. Biol. Chem.* **266**:18324.
14. Pinckard, J.K. *et al.* (1997) *J. Biol. Chem.* **272**:10784.
15. Engelmann, H. *et al.* (1990) *J. Biol. Chem.* **265**:1531.

MANUFACTURING SPECIFICATIONS

Animal-Free Manufacturing Conditions

Our dedicated controlled-access animal-free laboratories ensure that at no point in production are the products exposed to potential contamination by animal components or byproducts. Every stage of manufacturing is conducted in compliance with R&D Systems' stringent Standard Operating Procedures (SOPs). Production and purification procedures use equipment and media that are confirmed animal-free.

Production

- All molecular biology procedures use animal-free media and dedicated labware.
- Dedicated fermentors are utilized in committed animal-free areas.

Purification

- Protein purification columns are animal-free.
- Bulk proteins are filtered using animal-free filters.
- Purified proteins are stored in animal-free containers in a dedicated cold storage room.

Quality Assurance

- Low Endotoxin Level.
- No impairment of biological activity.
- High quality product obtained under stringent conditions.
- For ex vivo research or bioproduction, [additional documentation](#) can be provided.

[Please read our complete Animal-Free Statement](#)