

Lot #XXXXX

BostonBiochem

MATERIAL DATA SHEET

100X 20S Proteasome Activation Solution

Cat. # B-30

This buffer contains the optimal formulation for the chemical activation of all isoforms of constitutive 20S and immuno 20S proteasome enzymes including human (**E-360**), rabbit (**E-350**), rat (**E-375**), and mouse.

Product Information

Stock: 100 µl of 3% SDS

Use & Storage

Use: Dilute to 1X for a final concentration of 0.03% in assay. Dilute in reaction buffers prior to adding proteasome (do not add directly to solutions already containing proteasome). Use with peptide substrates containing arginine residues (**S-290**, **S-300**) is not recommended, as precipitation of the substrate may occur. In these instances non-chemical activation may provide better results.

Storage: Store tightly capped at room temperature. If stored at -20°C, make sure reagent is completely thawed and well-mixed before use.

For Laboratory Research Use Only, Not For Use in Humans

Rev: 12/07/2015

840 Memorial Drive, Cambridge, MA 02139 Phone: 617-576-2210 FAX: 617-492-3565
www.bostonbiochem.com

The contents of this datasheet (unless otherwise noted) are Copyright © 2008 Boston Biochem, Inc. All rights reserved. Duplication in whole or in part is strictly prohibited without the express written consent of Boston Biochem, Inc. "Boston Biochem" is a Trademark of Boston Biochem, Inc., 840 Memorial Drive, Cambridge, MA 02139.