



#### ORDERING INFORMATION

**Catalog Number:** BAF467

**Lot Number:** CSF02

**Size:** 50 µg

**Formulation:** 0.2 µm filtered solution in PBS with BSA

**Storage:** -20° C

**Reconstitution:** sterile 0.1% BSA in TBS

**Specificity:** mouse EphB2

**Immunogen:** NS0-derived rmEphB2 extracellular domain

**Ig Type:** mouse EphB2 extracellular domain specific goat IgG

**Application:** Western blot

## ***Biotinylated Anti-mouse EphB2 Antibody***

### ***Preparation***

Produced in goats immunized with purified, NS0-derived, recombinant mouse EphB2 (rmEphB2) extracellular domain. EphB2 specific IgG was purified by mouse EphB2 affinity chromatography and then biotinylated.

### ***Formulation***

Lyophilized from a 0.2 µm filtered solution in phosphate-buffered saline (PBS) containing 50 µg of bovine serum albumin (BSA) per 1 µg of antibody.

### ***Reconstitution***

Reconstitute with sterile Tris-buffered saline pH 7.3 (20 mM Trizma base, 150 mM NaCl) containing 0.1% BSA. If 1 mL of buffer is used, the antibody concentration will be 50 µg/mL.

### ***Storage***

Lyophilized samples are stable for twelve months from date of receipt when stored at -20° C to -70° C. Upon reconstitution, the antibody can be stored at 2° - 8° C for 1 month without detectable loss of activity. Reconstituted antibody can also be aliquotted and stored frozen at -20° C to -70° C **in a manual defrost freezer** for six months without detectable loss of activity. **Avoid repeated freeze-thaw cycles.**

### ***Specificity***

This antibody has been selected for use as a detection antibody in mouse EphB2 western blots. In western blots, this antibody shows approximately 5% cross-reactivity with rmEphB6, rrEphB1 and rmEphB3.

### ***Application***

**Western Blot** - This antibody can be used at 0.1 - 0.2 µg/mL with the appropriate secondary reagents to detect mouse EphB2. The detection limit for rmEphB2 is approximately 10 ng/lane under non-reducing and reducing conditions.

**Optimal dilutions should be determined by each laboratory for each application.**