

DESCRIPTION


Species Reactivity	Human/Mouse
Specificity	Detects human and mouse XIAP in Western blots.
Source	Polyclonal Goat IgG
Purification	Antigen Affinity-purified
Immunogen	<i>E. coli</i> -derived recombinant human XIAP Met1-Ser497 (Ser162Cys) Accession # P98170
Formulation	Lyophilized from a 0.2 µm filtered solution in PBS with BSA as a carrier protein. See Certificate of Analysis for details.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the *Technical Information* section on our website.

	Recommended Concentration	Sample
Western Blot	0.25 µg/mL	See Below

DATA



PREPARATION AND STORAGE

Reconstitution	Reconstitute at 0.2 mg/mL in sterile PBS.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	<p>Use a manual defrost freezer and avoid repeated freeze-thaw cycles.</p> <ul style="list-style-type: none"> ● 12 months from date of receipt, -20 to -70 °C as supplied. ● 1 month, 2 to 8 °C under sterile conditions after reconstitution. ● 6 months, -20 to -70 °C under sterile conditions after reconstitution.

BACKGROUND

XIAP (X-chromosome linked inhibitor of apoptosis) is a member of the apoptosis (IAP) family of proteins that inhibit caspases. The BIR2 domain of XIAP inhibits caspase-3 and caspase-7. The ability of XIAP to inhibit caspases is prevented by SMAC/Diablo through binding to XIAP-BIR2 and -BIR3 domains.