



***Affinity Purified
Donkey Anti-rabbit IgG (H & L) Antibody***

Catalog Number: D-301-C-ABS2

Lot Numbers: WBH0208051

WBH0208071

WBH0208091

Product Specifications

- | | |
|-----------------------|--|
| Amount | ◆ 500 µg |
| Form | ◆ Lyophilized |
| Purity | ◆ Affinity purified |
| Buffer | ◆ 10 mM sodium phosphate, 0.85% (w/v) sodium chloride, pH 7.2 (PBS) with 5% trehalose |
| Preservative | ◆ 0.05% (w/v) sodium azide |
| Reconstitution | ◆ Reconstitute with sterile PBS. If 0.5 mL of PBS is used, the antibody concentration will be 1 mg/mL. |
| Storage | ◆ 2° - 8° C |
| Specificity | ◆ This antibody reacts with heavy chains on rabbit IgG, and with light chains on all rabbit immunoglobulins. No reactivity observed with other serum proteins, or IgG from bovine, chicken, goat, guinea pig, horse, human, mouse, rat, sheep, or Syrian hamster, by IEP.
◆ This antibody adsorbed against solid phase human IgG, bovine IgG, mouse IgG, rat IgG, and goat IgG. |
| Note | ◆ Sodium azide may react with lead and copper plumbing to form explosive metal azides. Flush with large amounts of water during disposal. |