

*Affinity Purified*  
**Donkey Anti-sheep IgG (H & L) Antibody**

**Catalog Number:** D-501-C-ABS

**Lot Numbers:** WBM0319061

WBM0321031

***Product Specifications***

<b>Amount</b>	◆ 0.5 mg
<b>Form</b>	◆ Lyophilized
<b>Purity</b>	◆ Affinity purified
<b>Buffer</b>	◆ 10 mM sodium phosphate, 0.85% (w/v) sodium chloride, pH 7.2 (PBS) with 5% trehalose
<b>Preservative</b>	◆ 0.05% (w/v) sodium azide
<b>Reconstitution</b>	◆ Reconstitute with sterile PBS. If 0.5 mL of PBS is used, the antibody concentration will be 1 mg/mL.
<b>Storage</b>	◆ 2° - 8° C
<b>Specificity</b>	◆ This antibody reacts with heavy chains on sheep IgG, and with light chains on all sheep immunoglobulins. No reactivity is observed to non-immunoglobulin sheep serum proteins. This antibody is adsorbed against solid phase human, mouse, bovine, and rabbit serum and human and rabbit IgG.
<b>Note</b>	◆ Sodium azide may react with lead and copper plumbing to form explosive metal azides. Flush with large amounts of water during disposal.