

---

**MATERIAL DATA SHEET**

---

**Recombinant Human E6AP/UBE3A****Cat. # E3-230****Product Information**

<b>Quantity:</b>	50 µg
<b>MW:</b>	102 kDa
<b>Source:</b>	<i>Spodoptera frugiperda</i> , Sf21 (baculovirus)-derived Contains a C-terminal 6-His tag Accession # Q05086
<b>Stock:</b>	X mg/ml (X µM) in 50 mM HEPES pH 8.0, 100 mM NaCl, 1 mM TCEP
<b>Purity:</b>	>95%, by SDS-PAGE under reducing conditions and visualized by Colloidal Coomassie® Blue stain.

**Use & Storage**

<b>Storage:</b>	<b>Use a manual defrost freezer and avoid repeated freeze-thaw cycles.</b> <ul style="list-style-type: none"><li>• 12 months from date of receipt, -70 °C as supplied.</li><li>• 3 months, -70 °C under sterile conditions after opening.</li></ul>
-----------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

**Literature****References:**

*For research use only. Not for use in humans.*