

DESCRIPTION

Species Reactivity	Human
Specificity	Detects human BPIFB1 in sandwich ELISAs.
Source	Monoclonal Mouse IgG _{2B} Clone # 1023923
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	Human embryonic kidney cell HEK293-derived human BPIFB1 Thr22-Gln484 Accession # Q8TDL5
Conjugate	Alexa Fluor 594 Excitation Wavelength: 590 nm Emission Wavelength: 617 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide
*Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.	

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

ELISA	Optimal dilution of this antibody should be experimentally determined.
Immunohistochemistry	Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

BPIFB1 (BPI fold-containing family B member 1), also named LPLUNC1, is a member of the BPI-fold (BPIF) containing/Plunc (palate, lung, and nasal epithelium clone) superfamily of putative innate defense molecules which are predominantly expressed in regions of the oral cavity, nasopharynx and upper respiratory tract (1, 2). BPIF proteins exist as two subgroups, BPIFA (formally SPLUNCs) and BPIFB (formally LPLUNCs) (1, 3). BPIFA proteins have structural homology to the N-terminal domain of BPI whereas BPIFB proteins have structural homology to both domains of BPI (2). Human BPIFB1 cDNA encodes a 484 amino acid (aa) precursor protein with a putative 21 aa signal peptide and a 463 aa mature chain. The mature human BPIFB1 shares approximately 57% and 56% aa sequence identity with mouse and rat BPIFB1, respectively. BPIF proteins appear to exhibit distinct tissue and cell specific expression patterns with various family members, being localized to a number of glandular structures within the upper respiratory tract, nasopharyngeal regions and oral cavity where they are secreted from these tissues and are found in high levels in saliva and nasal and respiratory lining fluids (2). BPIFB1 plays a role in diverse functions, including neutralizing endotoxin (LPS) in septic shock patients, inhibition of endothelial cell growth, dendritic cell maturation, as an anti-angiogenic, chemoattractant or opsonization agent (2). Although less characterized as BPIFA1, BPIFB1 may also function as an innate immune molecule sensing and responding to Gram-negative bacteria (4). BPIFA1 and BPIFB1 expression was increased in late stage chronic obstructive pulmonary disease (COPD) patients, and elevated levels correlate with disease severity (5). BPIFB1 is also upregulated in cystic fibrosis (CF) lung disease and may play a role in the pathogenesis of the disease (6).

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.