

# Human VIGR/GPR126 Alexa Fluor® 594-conjugated Antibody

Monoclonal Mouse IgG<sub>2A</sub> Clone # 1044430

Catalog Number: FAB10577T

100 µg

## DESCRIPTION

|                           |                                                                                                                                                                                                                                                                                         |
|---------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Species Reactivity</b> | Human                                                                                                                                                                                                                                                                                   |
| <b>Specificity</b>        | Detects human VIGR/GPR126 protein in direct ELISAs.                                                                                                                                                                                                                                     |
| <b>Source</b>             | Monoclonal Mouse IgG <sub>2A</sub> Clone # 1044430                                                                                                                                                                                                                                      |
| <b>Purification</b>       | Protein A or G purified from hybridoma culture supernatant                                                                                                                                                                                                                              |
| <b>Immunogen</b>          | Human embryonic kidney cell HEK293-derived human VIGR/GPR126 protein<br>Cys38-Lys437<br>Accession # AAH75798.1                                                                                                                                                                          |
| <b>Conjugate</b>          | Alexa Fluor 594<br>Excitation Wavelength: 590 nm<br>Emission Wavelength: 617 nm                                                                                                                                                                                                         |
| <b>Formulation</b>        | Supplied 0.2 mg/mL in a saline solution containing BSA and Sodium Azide.<br><br>*Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions. |

## APPLICATIONS

**Please Note:** Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

|                       | <b>Recommended Concentration</b> | <b>Sample</b>                                             |
|-----------------------|----------------------------------|-----------------------------------------------------------|
| <b>Flow Cytometry</b> | 0.25-1 µg/10 <sup>6</sup> cells  | HEK293 Human Cell Line Transfected with Human VIGR/GPR126 |

## PREPARATION AND STORAGE

|                                |                                                                                                                                                    |
|--------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Shipping</b>                | The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.                                  |
| <b>Stability &amp; Storage</b> | <b>Protect from light. Do not freeze.</b> <ul style="list-style-type: none"> <li>12 months from date of receipt, 2 to 8 °C as supplied.</li> </ul> |

## BACKGROUND

VIGR (Vascular Inducible G Protein-coupled Receptor), also known as ADGRG6, DREG, and GPR126, is a neuronal 7 TM pass (G protein)-coupled receptor (GPCR) involved in myelination and glial and Schwann cell development (1, 2). Human VIGR cDNA encodes a 1221 amino acid (aa) residue membrane protein with a 37 aa signal peptide, a 825 aa extracellular domain (ECD) with 27 potential N-linked glycosylation sites, seven transmembrane segments that span between aa 863 and aa 1113, and a 108 aa residue cytoplasmic domain. Within ECD human VIGR shares 83% aa sequence identity with mouse and rat VIGR. VIGR is essential for the development of diverse organs (1, 2). Type IV collagen, a major constituent of the basement membrane, binds to VIGR and activates its signaling function (3). This interaction stimulated the production of cAMP in rodent Schwann cells, which require VIGR activity to differentiate, and in human embryonic kidney (HEK293) cells expressing exogenous VIGR. Laminin-211 binds a novel laminin-binding domain in VIGR N-terminal fragment between aa 446 and 807 (4). VIGR-Laminin-211 interactions regulate terminal differentiation and myelination by ensuring appropriate levels of cAMP for a given stage of Schwann cell development (4).

### References:

1. Ruggetti, A. *et al.* (2005) J. Immunol. **174**:7764.
2. Engelstaedter, V. *et al.* (2012) BMC Cancer **12**:600.
3. Taylor-Papadimitriou, J. *et al.* (1999) Biochim. Biophys. Acta **1455**:301.
4. Geng, Y. *et al.* (2012) Front Oncol. **2**:76.
5. Tanida, S. *et al.* (2013) J Biol Chem. **288**:31842.
6. Beatson, R. *et al.* (2016) Nat Immunol. **17**:1273.
7. Piyush, T. *et al.* (2017) Cell Death Differ. **24**:1937.

## PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc., and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.