

DESCRIPTION

Species Reactivity	HCoV-HKU1
Specificity	Detects Human coronavirus HCoV-HKU1 Spike RBD in direct ELISAs.
Source	Monoclonal Mouse IgG _{2B} Clone # 1045316
Purification	Protein A or G purified
Immunogen	Human embryonic kidney cell HEK293-derived HCoV-HKU1Spike RBD protein Accession # Q5MQD0.1
Conjugate	Alexa Fluor Plus 555 Excitation Wavelength: 558 nm Emission Wavelength: 572 nm
Formulation	Supplied 0.2 mg/mL in a saline solution containing BSA and Sodium Azide. *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

Western Blot Optimal dilution of this antibody should be experimentally determined.

DATA

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

HCoV-HKU1 was identified in Hong Kong in 2005 as a new human coronavirus (1). Coronaviruses are a family of viruses that are commonly comprised of a large plus-strand RNA genome and four structural proteins: Spike protein (S), Envelope protein (E), Membrane protein (M), and Nucleocapsid protein (N). There are two well-known human coronavirus families that infect humans: Alpha coronaviruses which includes HCoV-229E and HCoV-NL63; beta coronaviruses that includes HCoV-OC43, Severe Acute Respiratory Syndrome (SARS-CoV), Middle East Respirator Syndrome (MERS-CoV), and global pandemic Covid-19 (SARS-CoV2) (2). The HCoV-HKU1 Spike Protein (S Protein) is a glycoprotein that mediates membrane fusion and viral entry. As with most coronaviruses, proteolytic cleavage of the S protein generates two distinct peptides, S1 and S2 subunits. The S1 subunit is focused on attachment of the protein to the host receptor, while the S2 subunit is involved with cell fusion. The receptor binding domain (RBD) of HCoV-HKU1 is located at C-terminal region of S1 subunit, similar to SARS-COV, MERS-COV and SARS-COV2, but the RBD regions do not share significant amino acid sequence identity (3). HCoV-HKU1 has been demonstrated to bind specifically to 9-O-acetylated sialic acids (9-O-Ac-Sias) attached as terminal residues to glycan chains on glycoproteins and lipids, but additional receptors remain unknown (4). HCoV-HKU1, along with HCoV-OC43, differ from other coronaviruses in that their virions possess two types of surface projections, both involved in attachment: large "spikes" that are comprised of homotrimers of the S protein, and unique, smaller protrusions, comprised of the homodimeric hemagglutinin-esterase (HE) (5).

References:

1. Woo, P. *et al.* (2005) J. Virol. **79**:884.
2. Ogimi, C. *et al.* (2020) J Pediatric Infect Dis Soc doi: 10.1093/jpids/piaa037.
3. Qian, Z. *et al.* (2015) J. Virol. **89**:8816.
4. Huang, X. *et al.* (2015) J Virol **89**:7202.
5. Hulswit, R.J.G. *et al.* (2019) PNAS **116**:2681.

PRODUCT SPECIFIC NOTICES

This product is provided under an intellectual property license from Life Technologies Corporation. The transfer of this product is conditioned on the buyer using the purchased product solely in research conducted by the buyer, excluding contract research or any fee for service research, and the buyer must not (1) use this product or its components for (a) diagnostic, therapeutic or prophylactic purposes; (b) testing, analysis or screening services, or information in return for compensation on a per-test basis; or (c) manufacturing or quality assurance or quality control, and/or (2) sell or transfer this product or its components for resale, whether or not resold for use in research. For information on purchasing a license to this product for purposes other than as described above, contact Life Technologies Corporation, 5781 Van Allen Way, Carlsbad, CA 92008 USA or outlicensing@thermofisher.com.