bio-techne® RDSYSTEMS

DESCRIPTION	
Species Reactivity	Human
Specificity	Detects human FGFR1 alpha (IIIc) in direct ELISA.
Source	Monoclonal Mouse IgG _{2B} Clone # 1058809
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	Human embryonic kidney cell, HEK293-derived human FGFR1 alpha Arg22-Glu376 Accession # P11362.3
Conjugate	Alexa Fluor 647 Excitation Wavelength: 650 nm Emission Wavelength: 668 nm
Formulation	Supplied 0.2 mg/mL in a saline solution containing BSA and Sodium Azide.
	*Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS	
Please Note: Optimal dilutions should be de	termined by each laboratory for each application. General Protocols are available in the Technical Information section on our website.
Flow Cytometry	Titration recommended for optimal concentration with starting range of 0.1-1 μg/1 million cells. Sample used for this experiment was HEK293 cells transfected with hFGFR1 and eGFP.
Flow Cytometry	Titration recommended for optimal concentration with starting range of 0.1-1 μg/1 million cells. Sample used for th experiment was HEK293 cells transfected with hFGFR1 and eGFP.

PREPARATION AND STORAGE		
Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.	
Stability & Storage	Protect from light. Do not freeze.	
	 12 months from date of receipt, 2 to 8 °C as supplied. 	

BACKGROUND

Fibroblast growth factor receptor 1 (FGFR1) belongs to a family of type I transmembrane tyrosine kinases which mediate the biological functions of FGFs that are involved in a multitude of physiological and pathological cellular processes (1). The FGFR family is comprised of 4 structurally conserved members (FGFR1-4) all possessing an extracellular domain (ECD) with three immunoglobulin (Ig)-like domains, an acid-box region containing a run of acidic residues between the IgI and IgII domains, a transmembrane domain and cytoplasmic split tyrosine-kinase domain (1, 2). The ECD of mature, full-length FGFR1 shares 98% amino acid sequence identity with mouse FGFR1. Alternative splicing generates multiple forms of FGFR1-3, each with unique signaling characteristics (1-3). For FGFR1, alternative splicing of the ECD generates FGFR1A, FGFR1B, and FGFR1G isoforms of the receptor with the A isoform containing three Ig domains, while the B and G isoforms lack the N-terminal IgI domain (3). Additional splicing of the IgIII domain, results in IIIa, IIIb, or IIIc isoforms (3). Only the alpha isoform has been identified for FGFR3 and FGFR4 and FGFR4 also lacks the IIIb and IIIc splicing events (4). The FGFR splice variants also exhibit distinct and varying binding affinities for different FGF ligands (2). FGFRs mediate the FGF signaling cascade which regulate developmental processes including cellular proliferation, differentiation, and migration, morphogenesis, and patterning (5). FGFRs transduce the signals through three dominant pathways including RAS/MAPK, PI3k/AKT, and PLCY (6). FGFR1 the most abundant FGFR and is widely expressed in many adult human tissues, but the splice variants display distinct tissue-specific differences with IIIc splice variants expressed in mesenchymal tissue (4, 7, 8). Mutations in FGFR1 or misregulation of FGFR1 mediated signaling is found in multiple diseases, with FGFR1A(IIIc) specifically upregulated, from breast and pancreatic cancer to Pfeiffer syndrome and osteoarthritis (4, 9-11). A solubl

References:

- 1. Ornitz, D.M. and Itoh, N. (2015) Wiley Interdiscip. Rev. Dev. Biol. 4:215.
- 2. Zhang, X. et al. (2006) J Biol. Chem. 281:15694.
- 3. Ferguson, H.R. et al. (2021) Signaling Cells 10:1201.
- 4. Holzmann, K. et al. (2012) J Nucleic. Acids 2012:950508.
- 5. Xie, Y. et al. (2020) Sig. Transduct. Target Ther. 5:181.
- 6. Mossahebi-Mohammadi, M. et al. (2020) Front Cell Dev. Biol. 18:79.
- 7. Hughes, S.E. (1997) J. Histochem. Cytochem. 45:1005.
- 8. Delezoide, A.L. et al. (1998) Mech. Dev. 77:19.
- 9. Yamashita-Sugahara, Y. et al. (2016) Sci. Rep. 6:35908.
- 10. Teven, C.M. et al. (2014) Genes Dis. 1:199.
- 11. Babina, I. and Turner, N. (2017) Nat. Rev Cancer 17:318.

Rev. 1/19/2023 Page 1 of 2



Global bio-techne.com info@bio-techne.com techsupport@bio-techne.com TEL +1 612 379 2956 USA TEL 800 343 7475 Canada TEL 855 668 8722 China TEL +86 (21) 52380373 Europe | Middle East | Africa TEL +44 (0)1235 529449



Human FGFR1 alpha (IIIc) Alexa Fluor® 647-conjugated Antibody Monoclonal Mouse IgG_{2B} Clone # 1058809 Catalog Number: FAB11336R

100 µg

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.

Rev. 1/19/2023 Page 2 of 2



Global bio-techne.com info@bio-techne.com techsupport@bio-techne.com TEL +1 612 379 2956 USA TEL 800 343 7475 **Canada** TEL 855 668 8722 **China** TEL +86 (21) 52380373 **Europe | Middle East | Africa** TEL +44 (0)1235 529449