

#### DESCRIPTION

<b>Species Reactivity</b>	Human
<b>Specificity</b>	Detects human TAPBPR in direct ELISA.
<b>Source</b>	Monoclonal Mouse IgG <sub>2A</sub> Clone # 1059329
<b>Purification</b>	Protein A or G purified from hybridoma culture supernatant
<b>Immunogen</b>	Chinese Hamster Ovary cell line, CHO-derived human TAPBPR Ala19-Arg404 Accession # Q9BX59.2
<b>Conjugate</b>	Alexa Fluor 405 Excitation Wavelength: 405 nm Emission Wavelength: 421 nm
<b>Formulation</b>	Supplied 0.2 mg/mL in a saline solution containing BSA and Sodium Azide.  *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

#### APPLICATIONS

**Please Note:** Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

**Flow Cytometry** Titration recommended for optimal concentration with starting range of 0.1-1 µg/1 million cells. Sample used for this experiment was HEK293 cells transfected with Human TAPBL and eGFP vs irrelevant.

#### PREPARATION AND STORAGE

**Shipping** The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.

**Stability & Storage** **Protect from light. Do not freeze.**

- 12 months from date of receipt, 2 to 8 °C as supplied.

#### BACKGROUND

TAP-binding protein-like (TAPBPL), also known as TAP binding protein-related (TAPBPR) and Tapasin-related protein (TAPASINR) is a transmembrane protein of the Immunoglobulin (Ig) superfamily (1, 2). TAPBPR was originally isolated as a homologue to TAPASIN but more recently was identified as a novel B7 family-related molecule since it shares sequence, structural, and functional similarities to B7 family members (3). Mature human TAPBPR consists of a luminal domain containing an IgV and IgC domain, a transmembrane domain, and a cytoplasmic tail which lacks an ER retention motif. Within the luminal domain, mature human TAPBPR shares 70% and 71% amino acid sequence identity with mouse and rat TAPBPR, respectively. Multiple alternatively spliced TAPBPR isoforms are known to exist with unique properties (4). TAPBPR is widely expressed and, similar to TAPASIN, functions as both a chaperone protein and peptide editor of MHC class I, but in a peptide-loading complex (PLC) independent manner (5, 6). TAPBPR decreases the rate at which MHC class I molecules mature through the secretory pathway, a role which could be important for peptide selection by MHC class I molecules (7). TAPBPR is also expressed on the surface of T cells and antigen-presenting cells (APCs) and plays an inhibitory role in the proliferation and activation of T cells (4). TAPBPR can be expressed on various cancer cells like leukemia and has the potential to be used in the treatment of autoimmune diseases and transplant rejection, as well as cancer (4).

#### References:

1. Hermann, C. *et al.* (2015) *Tissue antigens* **85(3)**:155.
2. Teng, M. *et al.* (2002) *European Journal of Immunology* **32**:1059.
3. Lin, Y. *et al.* (2021). *EMBO Mol Med.* **13(5)**:13404.
4. Porter, K.M. *et al.* (2014) *Immunology* **142**:289.
5. Margulies, D. *et al.* (2020) *Current Opinion in Immunology* **64**:71.
6. Boyle, L.H. *et al.* (2013) *PNAS* **110**:3465.
7. Hermann, C. *et al.* (2013) *Journal of Immunology* **191(11)**:5743.

#### PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.