

## DESCRIPTION

|                           |                                                                                 |
|---------------------------|---------------------------------------------------------------------------------|
| <b>Species Reactivity</b> | Human                                                                           |
| <b>Specificity</b>        | Detects human GluR4 in direct ELISA.                                            |
| <b>Source</b>             | Monoclonal Mouse IgG <sub>2B</sub> Clone # 1059705                              |
| <b>Purification</b>       | Protein A or G purified from cell culture supernatant                           |
| <b>Immunogen</b>          | Synthetic peptide<br>Accession # P48058                                         |
| <b>Conjugate</b>          | Alexa Fluor 594<br>Excitation Wavelength: 590 nm<br>Emission Wavelength: 617 nm |
| <b>Formulation</b>        | Supplied 0.2 mg/mL in a saline solution containing BSA and Sodium Azide.        |

\*Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

## APPLICATIONS

**Please Note:** Optimal dilutions should be determined by each laboratory for each application. General Protocols are available in the Technical Information section on our website.

|                       |                                                                                                                                                              |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Flow Cytometry</b> | Titration recommended for optimal concentration with starting range of 0.1-1 µg/1 million cells. Sample used for this experiment was PBMC lymphocytes cells. |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|

## PREPARATION AND STORAGE

|                                |                                                                                                                   |
|--------------------------------|-------------------------------------------------------------------------------------------------------------------|
| <b>Shipping</b>                | The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below. |
| <b>Stability &amp; Storage</b> | <b>Protect from light. Do not freeze.</b><br>• 12 months from date of receipt, 2 to 8 °C as supplied.             |

## BACKGROUND

Metabotropic glutamate receptor 4 (mGluR4) is a 95-110 kDa Group III mGluR member of the Class C, GPCR 3 family of proteins. The mGluR family of receptors are 7-transmembrane (TM), non-ion channel glutamate-binding molecules that are linked to adenylyclase. mGluR4 is widely expressed, being found on/in multiple cell types, including colonic, breast, urinary and skin epithelium. It is also found on neurons where it is located presynaptically and regulates neurotransmitter release. In general, mGluR4 activation results in a slowing of either glutamine or GABA release. However, in discrete areas, it actually facilitates glutamine release. Mature human mGluR4 is an 880 amino acid (aa) 7-TM molecule (aa 33-912). It contains a 555 aa N-terminal extracellular region (aa 33-587), plus a 65 aa C-terminal cytoplasmic domain. mGluR4 potentially forms homodimers, and is reported to heterodimerize with both mGluR2 and mGluR8. Based on rodent, there is one 68 kDa isoform that is found in taste buds and shows an alternative start site at Met309. Two additional potential splice variants contain a second alternative start site at Met209, and an Asp substitution for aa 343-390, respectively. It is uncertain if human possesses a mGluR4b found in rodent. Over aa 1-518, human mGluR4 shares 96% aa identity with mouse mGluR4.

## PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.