

DESCRIPTION

Species Reactivity	Human
Specificity	Detects a synthetic peptide around aa 235 in Direct ELISA
Source	Monoclonal Mouse IgG _{2B} Clone # 1082626
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	ATP5A1 containing synthetic peptide Accession # P25705
Conjugate	Alexa Fluor 350 Excitation Wavelength: 346 nm Emission Wavelength: 442 nm
Formulation	Supplied 0.2mg/ml in 1X PBS with RDF1 and 0.09% Sodium Azide
*Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.	

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

Immunocytochemistry	Optimal dilution of this antibody should be experimentally determined.
Immunohistochemistry	Optimal dilution of this antibody should be experimentally determined.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied

BACKGROUND

ATP5F1A is a 553aa, 60kDa protein which encodes subunit alpha of mitochondrial ATP synthase (complex V), the rotary ATPase responsible for more than 90% of cellular ATP synthesis. F-type ATPases have 2 components, CF1 - the catalytic core - and CF0 - the membrane proton channel. CF1 has five subunits: alpha3, beta3, gamma1, delta1, epsilon1. CF0 has three main subunits: a, b and c. ATP5A1 is generally regarded to mediate tumor progression through mitochondrial signaling. It is overrepresented in an oxidative phosphorylation pathway associated with tumorigenesis and tumor progression. Dysregulation of this ATP synthase subunit is observed in human carcinomas.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.